

Table 2. Companies and hybrids included in the 2006 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables		
		Co.	MN	GRM	SRM				Trt	Co.	MN	GRM		SRM	Trt
AgriGold A6205Bt	CB-G	96	96	96	55	11	* Carharts Blue Top CR1010R	CB,RR-GJ	110	108	107	2	7*,12*,21		
* AgriGold A6225BtRR	CB,RR-GJ	98	98	100	55	6,11*	* Carharts Blue Top CR1857R	CB,CR,RR-GNJ	90	88	94	57	15*,17,22		
AgriGold A6275RWRR	CR,RR-NJ	100	100	97	55	11	* Carharts Blue Top CR1960R	CB,CR,RR-GNJ	100	100	98	57	6*,11,13*,16*,2		
AgriGold A6285BtRWRR	CB,CR,RR-GNJ	100	100	100	55	11	* Carharts Blue Top CR6585R	RR-J	85	87	92	2	15,17,22*		
AgriGold A6310Bt	CB-G	103	103	104	55	6	* Carharts Blue Top CR840RB	CB,RR-GJ	85	85	87	92	1	15,17*	
* AgriGold A6325RWRR	CR,RR-NJ	104	104	105	55	6*,12*	* Carharts Blue Top CX1956B	CB-G	95	95	93	1	11,16*		
AgriGold A6333BtRWRR	CB,CR,RR-GNJ	104	104	105	55	6,12	* Carharts Blue Top CX1958P	CB,CR-GN	95	96		57	11,13*,16		
AgriGold A6394Bt	CB-G	107	107	108	55	7,12									
AgriGold A6395BtRWRR	CB,CR,RR-GNJ	107	107	107	55	7,12									
* AgriGold XA6501BtRR	CB,RR-GJ	97	97		55	11*	* Cornelius C327YG	CB-G	103	103		55	6*		
							* Cornelius C382YG	CB-G	104	102	104	55	6,18*		
* Asgrow RX668RR2YGCB	CB,RR-GJ	107	110	108	50	7*	Cornelius C447RR	RR-J	106	108		55	7		
Asgrow RX715RR2YGCB	CB,RR-GJ	111	111		50	7	* Cornelius C462		106	106		55	7*		
							Cornelius C466Hx	CB,LL-LD	105	105		55	6		
* Blue River 26K21		88	88		87	96	22*	Cornelius C542RR	RR-J	108	108	110	55	7,19	
Blue River 34C17		94	93	94		3	10	* Cornelius C595		109	108	109	57	7*,19	
* Blue River 42A32		98	97	96		3	10*	Cornelius C626YGPL	CB,CR-GN	112	109		55	7	
Blue River 49M37		104	102		104	96	21	Cornelius C635YG	CB-G	112	110		55	7	
								Cornelius C684RRYG	CB,RR-GJ	112	110	110	55	7,19	
Brown 1990		85	84	85		108	17	Croplan Genetics 229RRBt	CB,RR-GJ	84	84	86		98	15,17
Brown 2830		85		87		108	17	Croplan Genetics 238TS	CB,CR,RR-GNJ	84	85	87		37	15,17
* Brown 3000YGCB	CB-G	90	90	88	92	48	15*,17	Croplan Genetics 253RRBt	CB,RR-GJ	87	87	88		37	15
Brown 3006		90	90	89		108	15	* Croplan Genetics 294RRBt	CB,RR-GJ	90	90		92	37	22*
Brown 4110BtLL	CB,LL-KD	92	95	95		109	16	* Croplan Genetics 296TS	CB,CR,RR-GNJ	90	90	88		37	15*
* Brown 4250YGCB	CB-G	95	95	94		43	16*	* Croplan Genetics 355TS	CB,CR,RR-GNJ	95	95	97		37	11,13*,14,16
Brown 4989RRLfy	Leafy,RR-J	96	98		100	110	23	* Croplan Genetics 364CRWRR	CR,RR-NJ	96	96	94	92	37	16*,
Brown 5232RRYGPlus	CB,CR,RR-GNJ	100	101	100		97	6,11	* Croplan Genetics 388RRBt	CB,RR-GJ	99	99	98		37	11*,14*,16*
* Brown 5310Hx1LL	CB,LL-LD	102	102	102		112	6*	* Croplan Genetics 421RRBt	CB,RR-GJ	100	100	100		37	6*,9*,11*,14*,16*
* Brown 5636		100	102		99	110	23*	* Croplan Genetics 491TS	CB,CR,RR-GNJ	103	100	102		37	6*,8,9*,12*,13*,14
Brown 5636RR2YGCB	CB,RR-GJ	102	101	103		111	6	* Croplan Genetics 566TS	CB,CR,RR-GNJ	105	106	104	105	37	6*,8*,9*,12*,13*,14*,16
* Brown 6810BtLL	CB,LL-KD	108		108		113	7,12*	* Croplan Genetics 591CRWRR	CR,RR-NJ	108	109	108	107	37	7*,8*,9*,19*
								Croplan Genetics 631TS	CB,CR,RR-GNJ	110			108	37	21
Brunner EXP102Bt	CB-G	102		102		54	12	* Croplan Genetics 645TS	CB,CR,RR-GNJ	111	110			37	7*,8*,9
Brunner EXP105Bt	CB-G	105			105	54	21	* Croplan Genetics 691BtLL	CB,LL-KD	112	112	109	109	37	7,8*,9,19,21
* Brunner EXP106		106			107	54	21*	Croplan Genetics 721TS	CB,CR,RR-GNJ	114			110	37	19,21
* Brunner EXP110		110		110		54	7*	* Croplan Genetics DS107RR2	RR-J	107	107		108	37	19*,21,
* Brunner EXP94	RR-J	94		95		54	16*								
Brunner EXP9506		95		95		3	10	* Crows 2121S	CB,RR-GJ	101	101	103	105	37	6*,12*,14*,2
Brunner EXP97Bt	CB-G	97		96	98	54	16,23	* Crows 4222S	CR,RR-NJ	107	107	108	109	37	7,9,19*
* Brunner OR1005		100	100	98		3	10*	* Crows 4843X	CR,RR-NJ	110		110	111	37	7,8*,9,19*
* Brunner OR1053		105	105	102		3	10*								
Brunner OR9004		90	90	94		3	10	Dahlco 2146		84		93		3	10
Brunner OR9406		95		97		3	10	Dahlco 2147	CB,RR-GJ	84	83	87		54	15
Brunner S1137Bt	CB-G	113		111		54	7	Dahlco 2395		92	93	95		3	10
Brunner S2065RR	RR-J	83		87		54	17	* Dahlco 2482	CB,RR-GJ	95	96	95		54	16*
Brunner S3403RRBt	CB,RR-GJ	84	84	87		54	17	Dahlco 2482		96	96	97		3	10
Brunner S3503CRWBt	CB,CR-GN	92	92	96		54	16	* Dahlco 2505	CB,RR-GJ	103	103	103		54	12*
* Brunner S3704Bt	CB-G	97	97	94		54	16*	Dahlco 3920	CB,RR-GJ	92	92	95		54	16
* Brunner S3792RRBt	CB,RR-GH	93	93	86		54	17*	Dahlco 3930	CB,RR-GJ	93	93		54	16	
Brunner S5086Bt	CB-G	102		103		54	12								

† Trait-Gene=Code: bmr-bm3=B; IMI:IT=C; LL-T25=D; Bt-CB-Bt176=F,Mon810=G,Bt11=K,TC1507=L; RR-MonGA21=H,Nk603=J,SytGA21=Q; Leafy-lfy=M; Bt-CR-Mon863=N, DAS591227=P; Unknown-X.

Table 2. Companies and hybrids included in the 2006 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	
		Co.	MN	GRM	SRM				Trt	Co.	MN	GRM		SRM
Dahlco 4013	CR,RR-NJ	101	101	54		12	* Dekalb DKC48-53(RR2YGCB)	CB,RR-GJ	98	100	96	99	7	11,14,16
Dahlco 4851Hx	CB,LL-KD	85	88	54		15	Dekalb DKC50-20(RR2YGCB)	CB,RR-GJ	100	100	104	98	55	20
* Dahlco X5003		100	98	3		10*	* Dekalb DKC51-39(RR2YGPL)	CB,CR,RR-GNJ	101	100	102	100	7	6,8*,9*,12
* Dahlco X6003	CR-N	100	102	54		11*	* Dekalb DKC52-40(RR2YGPL)	CB,CR,RR-GNJ	102	100	101		7	6*,8*,9*
							* Dekalb DKC52-63(RR2YGCB)	CB,RR-GJ	102	100	102		7	6*,12*
Dahlman D4215CB	CB-G	85	85	87	54	15	* Dekalb DKC54-53(RR2YGCB)	CB,RR-GJ	104	105	106	106	7	18*,21
* Dahlman D4515	CB-G	92	92	91	54	16*	Dekalb DKC55-82(RR2)	RR-J	105	105		106	7	21
* Dahlman D4530		90	86	54		15*	* Dekalb DKC57-79(RR2YGPL)	CB,CR,RR-GNJ	107	105	108		7	7*,8*,9*
* Dahlman D4801		96	96	93	54	16*	Dekalb DKC58-19(RR2)	RR-J	108	109			7	7
* Dahlman D4815	CB-G	95	95	95	54	16*	* Dekalb DKC58-78(YGCB)	CB-G	108	108	106		54	7,8*,9
Dahlman R4215	RR-J	85	85	86	54	15	Dekalb DKC58-80(RR2YGCB)	CB,RR-GJ	108	108		105	7	21
Dahlman R4215CB	CB,RR-GJ	85	85	87	54	15	* Dekalb DKC60-18(RR2YGPL)	CB,CR,RR-GNJ	110	108			7	7*,8*,9*
* Dahlman R4515	RR-J	90	90	87	54	15*								
* Dahlman R4515YGP	CB,CR,RR-GNJ	92	92	95	99	16*	Dyna Gro 50K07	RR-J	78	78	83		100	17
* Dahlman R45S15	CB,RR-GJ	92	92	94	54	16*	Dyna Gro 51P33	CB,RR-GJ	85	84	86		100	17
* Dahlman R4807CB	CB,RR-GJ	96	96	96	54	16*	* Dyna Gro 53F09	CB-G	93	95	95		100	11,16*
* Dahlman R4815	RR-J	95	95	94	54	16*	* Dyna Gro 53P30	CB,RR-GJ	92	95	92	97	100	16,17,2
							Dyna Gro 55B02	CB,CR,RR-GNJ	103	102	102		100	12
Dairyland DST10563		100		101	33	23	Dyna Gro 55B65	CB,CR,RR-GNJ	105		104		100	6
Dairyland HiDF3002		102		106	33	21	* Dyna Gro 55F43	CB-G	105	104	104	100	100	6*,12*,21*
* Dairyland HiDF3007		109		109	33	19*,21	* Dyna Gro 55P41	CB,RR-GJ	101	102	99		100	16*
Dairyland HiDF3087RR	RR-J	87		92	33	24	* Dyna Gro 56P22	CB,RR-GJ	110	108	106	100		7*,21
Dairyland HiDF3094		94		98	33	23								
Dairyland Stealth 1600		100	99		33	6,11	* Farm Advantage 5406	CB,LL-LD	105	106	104		43	6*
* Dairyland Stealth 1603		103	102		33	6*	Farm Advantage 9699L	CB-K	99	99	98		43	11
* Dairyland Stealth 1612		110		111	33	19*								
Dairyland Stealth 1615		113		112	33	19	Foundation Direct 8797O		97	93		3		10
Dairyland Stealth 1806		106	108		33	7								
* Dairyland Stealth 4006	CR,RR-NJ	106	108		33	7*	* Garst 8533YPL	CB,CR-GN	108	108			37	7*
* Dairyland Stealth 4009	CR,RR-NJ	111	110		33	7*	Garst 8534YG1GT	CB,RR-GP	101	108		106	37	21
* Dairyland Stealth 5007	CB-G	107		110	33	19*	* Garst 8590GT	RR-P	100		100		37	20*,23
* Dairyland Stealth 5194	CB-G	94	94	95	33	16*	* Garst 8677GT	RR-P	105	103	106		37	6*,21
Dairyland Stealth 5201	CB-G	101	101	102	33	6,12	* Garst 8688GT	RR-P	104	102	105		37	6,21*
* Dairyland Stealth 5204	CB-G	104	104	103	33	6*,12*	Garst 8827CBLLGT	CB,LL,RR-KP	96	98			37	11
* Dairyland Stealth 5497	CB-G	98	97	95	33	11,16*	* Garst 8866RR	RR-J	90		94		37	22*
* Dairyland Stealth 5503	CB-G	105	106	104	33	6,12*	* Garst 8880YG1	CB-G	96	96	96		37	11,16*
Dairyland Stealth 5611	CB-G	112	110		33	7	* Garst 8921YG1RR	CB,RR-GJ	92	90	91		37	16*,17
* Dairyland Stealth 7184	CB,RR-GJ	85	87		33	15,17*	Garst 8980YG1	CB-G	84	84	87		37	15,17
Dairyland Stealth 7191	CB,RR-GJ	91	91	90	33	16,17								
* Dairyland Stealth 7196	CB,RR-GJ	96	96		33	11,16*	Gold Country 8301CBR	CB,RR-GJ	83	83	85		55	15
Dairyland Stealth 7201	CB,RR-GJ	101	103		33	12	Gold Country 8902R	RR-J	89	89	87		55	15
Dairyland Stealth 7685	CB,RR-GJ	86	87		33	15,17	Gold Country 9606SLS		96	96		98	55	23
Dekalb DKC33-11(RR2YGCB)	CB,RR-GJ	83	85	87	7	17	Golden Harvest EX68636Hx	CB-L	109	109		102		7
Dekalb DKC33-72(RR2)	RR-J	83	85	83	7	17	* Golden Harvest H6466CBGT	CB,LL,RR-KP	88	86	87	102		15,22*
Dekalb DKC35-02(RR2YGCB)	CB,RR-GJ	85	85	86	7	15,17	* Golden Harvest H7625HxLL	CB,LL-L	99	101	100	102		6*,11*,16*
* Dekalb DKC40-08(RR2YGCB)	CB,RR-GJ	90	90	89	92	7	Golden Harvest H7935Hx	CB-L	103	103	103		102	6,12
* Dekalb DKC42-88(RR2YGPL)	CB,CR,RR-GNJ	92	90	98	7	23*	Golden Harvest H8265CB	CB,LL-K	107	105	107		102	7,12
* Dekalb DKC44-92(RR2)	RR-J	94	95	97	7	11,16*	* Golden Harvest H8445CB	CB,LL-K	106	107	108	109	102	7,19*
* Dekalb DKC45-82(RR2)	RR-J	95	95	95	7	11,16*	* Golden Harvest H8713CB	CB-K	111	109	109	108	102	7,19*

† Trait-Gene=Code: bmr-bm3=B; IMI:IT=C; LL-T25=D; Bt-CB-Bt176=F, Mon810=G, Bt11=K, TC1507=L; RR-MonGA21=H, Nk603=J, SytGA21=Q; Leafy-lfy=M; Bt-CR-Mon863=N, DAS591227=P; Unknown-X.

Table 2. Companies and hybrids included in the 2006 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	
		Co.	MN	GRM	SRM				Trt	Co.	MN	GRM		SRM
Golden Harvest H8959Bt	CB-G	111		110	102	19	* Jung 7437RRYGPL	CB,CR,RR-GNJ	96	94	97		95	6,11,16*
Golden Harvest H9205Hx	CB-L	112		110	102	19	Jung 7711RRYGPL	CB,CR,RR-GNJ	111	112	106		95	7
Great Lakes 3663Bt	CB-G	86	87	55		15	Kaltenberg K2405RRBt	CB,RR-GJ	81		87		93	17
* Great Lakes 4041		90	87	55		15*	* Kaltenberg K2801RR	RR-J	84	86	92	93		15,17
* Great Lakes 4297BtRR	CB,RR-GJ	92	93	94	55	16*	Kaltenberg K3535Bt	CB-G	89		88		93	17
* Great Lakes 4415BtRR	CB,RR-GJ	94	95	55		16*	Kaltenberg K3914Bt	CB-G	92	91	88		93	17
* Great Lakes 4689BtRR	CB,RR-GJ	96	99	55		11*	* Kaltenberg K3919RRBt	CB,RR-GJ	92	92	95	96	93	16*,23
Great Lakes 5098G3	CB,CR,RR-GNJ	100	99	55		11	* Kaltenberg K4012RRBt	CB,RR-GJ	97	97	96		93	11*,16*
Great Lakes 5377G3	CB,CR,RR-GNJ	103	104	55		12	* Kaltenberg K4353RRPlus	CB,CR,RR-GNJ	100		99		93	16*
Great Lakes 5416RWRR	CR,RR-NJ	104	105	104		12	* Kaltenberg K4666RR	RR-J	96	96	91		93	16*,17*
							* Kaltenberg K4688Bt	CB-G	96	96	95		93	16*
Growmark FS3327RRYGCB	CB,RR-GJ	93	95	103		16	Kaltenberg K4826PLUS	CB,CR-GN	97		94		93	11,16
* Growmark FS4366YGCB	CB-G	102	102	103	99	12*	Kaltenberg K5215Bt	CB-G	102	102	103	102	93	12,16,2
* Growmark FS4475YGCB	CB-G	104	104	55		12*	Kaltenberg K5539LLBtHx	CB,LL-L	105		106		93	6
Growmark FS5116NDRR	RR-J	105		107	105	21	Kaltenberg K6231Bt	CB-G	107		109		93	7
* Growmark FS5197RRYGPL	CB,CR,RR-GNJ	105	105	99		12*	* Kaltenberg K6573RRPlus	CB,CR,RR-GNJ	109		108		16	7*
Growmark FS5527RRYGCB	CB,RR-GJ	107	109	103		7	Kaltenberg K7221Bt	CB-G	112		109		93	7
Growmark FS5586Hx1	CB,LL-L	107	109	99		7	* Kaltenberg K8099LFRR	Leafy,RR-J	99			92	93	24*
* Growmark FS6297RRYGCB	CB,RR-GJ	110	110	110	99	7,19*	Kaltenberg K8105LFRR	Leafy,RR-J	105	105		106	93	23
* Growmark FS6317NDYGCB	CB-G	109		110	106	19*	Kaltenberg K8110LF	Leafy-	110		105		93	19,23
* Growmark FS6427RRYGCB	CB,RR-GJ	111	109	110	99	7*,19*	Kaltenberg K8112LF	Leafy-	112		108		93	21,23
Hughes 2946R+	CB,CR,RR-GNJ	101	101	55		6	Kruger K0109		109	108			7	7
* Hughes 3235CB	CB-G	101	101	55		6*	Kruger K0186		86	85		102		15
* Hughes 3371R+	CB,CR,RR-GNJ	103	102	55		6*	* Kruger K0191		91	94		102		16*
Hughes 3836CBL	CB,LL-GD	104	103	43		6	Kruger K0192		92	95		7		16
* Hughes 5401R+	CB,CR,RR-GNJ	107	108	43		7*	Kruger K0309		109	108		7		7
* Hughes 6226R+	CB,CR,RR-GNJ	109	109	55		7*	* Kruger K0404		104	103		7		6*,12
* Hughes 6951R+	CB,CR,RR-GNJ	110	111	55		7*	Kruger K0408		108	108		7		7
							Kruger K0409		109	108		102		7
Jeske JS636		95	95	1		16	Kruger K0508		108	105	106	102		7
* Johnson Seeds 2194RRBt	CB,RR-GJ	90	94	87	7	15*,17*	* Kruger K0510A		110	109		7		7*
* Johnson Seeds 4102Bt	CB-G	102	100	7		12*,16*	* Kruger K0510B		110	110		102		7*
* Johnson Seeds 4296RRBt	CB,RR-GJ	95	94	7		16*	Kruger K0599		99	98		102		11
Johnson Seeds 4651RRBt	CB,RR-GJ	102	102	93	7	12	* Kruger K0605B		105	105	103	102		6*,12*
* Johnson Seeds 4747RRBt	CB,RR-GJ	98	95	7		16*	Kruger K0610		110	108		7		7
* Johnson Seeds 4855Hx	CB-L	86	85	89	7	15,17*	Kruger K0692		92	92	93	102		16
Johnson Seeds 4897YGPLRR	CB,CR,RR-GNJ	102	102	102	7	12	Kruger K1190RR	RR-J	90	86		7		15
* Johnson Seeds 4986Bt	CB-G	98	97	7		16*	Kruger K1195RR	RR-J	95	95		7		11,14,16
* Johnson Seeds 4990Bt	CB-G	90	86	7		15*,17*	Kruger K1292RR	RR-J	92	95		7		16
* Jung 3432	CB-G	95	95	95	95	16*	* Kruger K1500RR	RR-J	100	100	98	102		6*,9*,11*,14*,16*
* Jung 4220	CB,RR-GJ	88	88	87	95	15*	Kruger K1503RR	RR-J	103	102		102		12
Jung 4440	CB,RR-GJ	95	95	97	95	11	Kruger K1584RR	RR-J	84	84	86	102		15
Jung 6572YGCB	CB-G	105	103	95		6	Kruger K1587RR	RR-J	87	86		7		15
Jung 7202	CB,CR,RR-GNJ	83	87	95		15	* Kruger K1606RR	RR-J	106	107		7		7*,9*,12*,14*
Jung 7418A	CB,CR,RR-GNJ	92	92	95	95	16	* Kruger K1700RR	RR-J	100	98		7		6,11,16*
Jung 7422RRYGPL	CB,CR,RR-GNJ	98	98	100	95	6,11	Kruger K2105RRYGCB	CB,RR-GJ	105	105		102		6
							* Kruger K2288RRYGCB	CB,RR-GJ	88	87		7		15*
							Kruger K2301RRYGCB	CB,RR-GJ	101	102		102		12

† Trait-Gene=Code: bmr-bm3=B; IMI:IT=C; LL-T25=D; Bt-CB-Bt176=F,Mon810=G,Bt11=K,TC1507=L; RR-MonGA21=H,Nk603=J, SvtGA21=Q; Leafy-lfy=M; Bt-CR-Mon863=N, DAS591227=P; Unknown-X.

Table 2. Companies and hybrids included in the 2006 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	
		Co.	MN	GRM	SRM				Trt	Co.	MN	GRM		SRM
* Kruger K2499RRYGCB	CB,RR-GJ	99	100	100	7	6*,9,11*,14*,16*	LG Seeds LG2517Bt	CB-G	103	103	106	62	12	
* Kruger K2506RRYGCB	CB,RR-GJ	106	106	107	102	7,9*,12,14	LG Seeds LG2526SP		105	105	107	7	21	
* Kruger K2509RRYGCB	CB,RR-GJ	107	106	102	102	7,9*,12*,14*	LG Seeds LG2533		105	110	104	7	6	
Kruger K2511RRYGCB	CB,RR-GJ	111	110		7	7	LG Seeds LG2540Bt	CB-G	109	110	108	7	7	
Kruger K2585RRYGCB	CB,RR-GJ	85	85	87	102	15	LG Seeds LG2545Bt	CB-G	107	110		7	7	
Kruger K2688RRYGCB	CB,RR-GJ	88	88	87	102	15	* LG Seeds LG2553		109		107	7	19*	
* Kruger K2697RRYGCB	CB,RR-GJ	97	97	96	102	11*,14,16*	* LG Seeds LG2633Bt	CB-G	113		110	7	19*	
Kruger K5209BtLL	CB,LL-K	109	109		102	7								
Kruger K5310YGCB	CB-G	110	110		102	7	* Laser L7H07Bt	CB-G	96	96	95	97	102	11,16*,20
Kruger K5406BtLL	CB,LL-K	106	110		102	7	* Laser L8H69Bt	CB-G	103		99	102		23*
* Kruger K5407BtLL	CB,LL-K	107	109		102	7*								
Kruger K5497YGCB	CB-G	97	96		7	11	Legacy Seeds L2020RRBt	CB,RR-GJ	82		87		57	15
* Kruger K5498YGCB	CB-G	98	98		102	11*	* Legacy Seeds L2927RRBt	CB,RR-GJ	92	96	95	57		16,23*
* Kruger K5504YGCB	CB-G	103	104	103	102	6*,12*	* Legacy Seeds L3077RRCRW	CB,CR,RR-GNJ	95	96	98	57		11,23*
Kruger K5505YGCB	CB-G	105	105	104	102	6,12	* Legacy Seeds L3285RRBt	CB,RR-GJ	96	97	99	57		11,16,20
* Kruger K5588BtLL	CB,LL-K	88	85		102	15*	* Legacy Seeds L3295RRBt	CB,RR-GJ	97	96		57		11*
* Kruger K5589BtLL	CB,LL-K	89	87		102	15*	* Legacy Seeds L4267RRBt	CB,RR-GJ	102		103	57		18*,21
* Kruger K5593BtLL	CB,LL-K	93	94		102	16*	Legacy Seeds L4660RRBt	CB,RR-GJ	102	103		57		6,12
* Kruger K5596	CB,LL-K	96	95		102	11*	Legacy Seeds L4840RRBt	CB,RR-GJ	103	101		57		6
Kruger K5597BtLL	CB,LL-K	97	96		102	11	* Legacy Seeds L6000BMR	BMR-B	110		112	57		19*
* Kruger K6503TS	CB,CR,RR-GNJ	103	103	102	102	6*,8,12,13*								
Kruger K6585TS	CB,CR,RR-GNJ	85	85	87	102	15	* Lemke 3081Bt	CB-G	90		87		56	15*
* Kruger K6603TS	CB,CR,RR-GNJ	103	104		7	6*,8,12,13*	* Lemke 3087RRPlus	CB,CR,RR-GNJ	90		88		57	15*
Kruger K6607TS	CB,CR,RR-GNJ	107	105		7	7,8,13	* Lemke 4011Bt	CB-G	95		95		2	16*
* Kruger K6611TS	CB,CR,RR-GNJ	111	109		7	7*	* Lemke 4075Plus	CB,CR-GN	97		96		57	16*
Kruger K9111YGCB	CB-G	111	111	108	102	7	Lemke 5087RRPlus	CB,CR,RR-GNJ	100	100		57		11
* Kruger K9203RRYGCB	CB,RR-GJ	103	102	104	102	6,9*,12,14	* Lemke 6061Bt	CB-G	105	102		2		12*
* Kruger K9212TS	CB,CR,RR-GNJ	112	109		102	7,8*								
Kruger K9310TS	CB,CR,RR-GNJ	110	109		102	7,8	* Midwest Genetics 4S502	CB,RR-GJ	97	97	99	98	37	11*,14*,20*
* Kruger K9392RRYGCB	CB,RR-GJ	92	93	94	102	11,16*	* Midwest Genetics 70503S	CB,RR-GJ	101	101	103	103	37	6*,9*,12*,14*,18*
* Kruger K9392TS	CB,CR,RR-GNJ	92	92	96	102	11,13*,16*	* Midwest Genetics 76172S	CB,RR-GJ	107	107	108	107	37	7,9,19*
* Kruger K9407TS	CB,CR,RR-GNJ	107	107	106	102	7,8,13*	* Midwest Genetics 77124X	CR,RR-NJ	110	109	109	37		7*,9,19*
* Kruger K9496		96	96	94	102	16*	* Midwest Genetics 7S322	CB,RR-GJ	111		110	37		19*
* Kruger K9496RR	RR-J	96	96	96	102	11,14,16*								
* Kruger K9496YGCB	CB-G	96	96	96	102	11,16*	* Mycogen 2A517	CB,LL-LD	101	101	102	107		6*,12*
Kruger K9593RRYGCB	CB,RR-GJ	93	95	94	102	11,16	* Mycogen 2D322	RR-J	92	91	94	107		16*
							Mycogen 2D326	CB,CR,RR-GNJ	92	92	95	107		16
* Kussmaul SB2983RRYGPlus	CB,CR,RR-GNJ	83	87		95	17*	Mycogen 2D468	CB,CR,RR-GNJ	99	99	97	107		11
Kussmaul SB414RR	RR-J	114	110		95	7	* Mycogen 2D555	CB,LL-LD	103	103	103	107		6,12*
Kussmaul SB500RR	RR-J	100	101	99	95	11	Mycogen 2D673	CB,LL-LD	109	109	108	107		7
* Kussmaul SB5595RR	RR-J	95	95	95	95	16*	Mycogen 2G677	CB,LL-LD	109	110		107		7
Kussmaul SB590RRYGPlus	CB,CR,RR-GNJ	90	89		95	17	Mycogen 2J271	RR-J	87	87	86	107		15,17
* Kussmaul SB597RRYGPlus	CB,CR,RR-GNJ	97	95		95	16*	* Mycogen 2J527	CB,LL-LD	103	103	104	107		6*,12
Kussmaul SB602RRYGPlus	CB,CR,RR-GNJ	102	103		95	12	Mycogen 2K546	CB,CR,RR-GNJ	103	103	102	107		6,12
* Kussmaul SB609RR	RR-J	109	107		95	7*	* Mycogen 2P174	CB,RR-GJ	85	85	87	107		15,17*
Kussmaul SB612RR	RR-J	112	112		95	7	* Mycogen 2R421	CB,RR-GJ	95	96	95	107		16*
							* Mycogen 2R426	CB-G	96	96	95	107		11,16*
* LG Seeds LG2475BtRR	CB,RR-GJ	96	97	96	7	16*	Mycogen 2R570	CB,LL-LD	104	104	103	107		6
* LG Seeds LG2496BtRR	CB,RR-GJ	100	100	99	7	11*	* Mycogen F2F444	BMR-B	102	98		107	107	21*
LG Seeds LG2507RR	RR-J	101	101	102	7	12	* Mycogen F2F485	BMR-B	98		98	107		20*,23

† Trait-Genes=Code: bmr-bm3=B; IMI:IT=C; LL-T25=D; Bt-CB-Bt176=F,Mon810=G,Bt11=K,TC1507=L; RR-MonGA21=H,Nk603=J, SvtGA21=Q; Leafy-lfy=M; Bt-CR-Mon863=N, DAS591227=P; Unknown-X.

Table 2. Companies and hybrids included in the 2006 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables		
		Co.	MN	GRM	SRM				Trt	Co.	MN	GRM		SRM	Trt
* Mycogen F2F566	BMR-B	105	82	105	107	18*,21*	OBrien OB1005		100	102		55	6		
* Mycogen F2F631	BMR-B	110	82	111	107	19*,21*	OBrien OB1077		107	108		55	7		
* Mycogen F2F699	BMR,CB,LL-BLD	112	85	113	107	19*	* OBrien OB1107		110	110		55	7*		
Mycogen TMF2H308	CB,LL-LD	96	96	99	107	20	* OBrien OB1113Bt	CB-G	110	109	107	55	7*,19*		
* Mycogen TMF2N422	RR-J	94		100	107	20*,23	OBrien OB955		95	95		55	11		
* Mycogen TMF2Q714	CB,LL-LD	110		109	107	19*,21	* OBrien OB966		96	97	97	55	11,20*		
* Mycogen TMF2Q731	CB,LL-LD	112		112	107	19*	OBrien OB977	CB-K	97	97		55	11		
* Mycogen X25102	CB,LL-LD	83	87		107	15,17*	OBrien OBX107		107	107		55	7,12		
* NK Brand N16-M1	CB,LL-KD	82	82	86	37	17*	* Organic UTC		103	102		3	10*		
* NK Brand N17-R3	CB,LL-KD	82	82	86	37	17*	* Organic UTC-Hand Weed		103	102	101	3	10*		
* NK Brand N20-R7	CB,LL-KD	85	85	86	37	17*	Partners Brand 479RRYGCB	CB,RR-GJ	100	104		95	6		
NK Brand N21-L4		86	86	86	37	17	Partners Brand 520Hx1	CB,LL-L	101	103		95	6		
* NK Brand N22-C2	CB,LL,RR-KDH	87	87	85	37	17*	* Partners Brand 578		107	110		95	7*		
* NK Brand N23-F7	CB,LL-KD	88	88	86	37	15,17*	Partners Brand 610YGCB	CB-G	110	109		95	7		
* NK Brand N27-W8	CB,LL-KD	91	91	90	92	37	16,17*								
NK Brand N29-A2	CB,LL-KD	92	92	94	37	16									
NK Brand N31-P2	Leafy,LL-D	92	102		92	37	24	* Pioneer 33N29		113	111	110	50	19*	
NK Brand N33-H6	Leafy-	93	101		97	37	23,2	* Pioneer 33T56		108	112	109	50	19*	
NK Brand N34-F1	CB,LL-KD	94	94	97	37	11	11	* Pioneer 34A18	CB,CR,LL-LP	109	109	109	106	50	7,19*,21*
NK Brand N35-B8	CB,LL-KD	96	96	98	37	6,11	6,11	* Pioneer 34A86		106	108		103	50	21*,23
* NK Brand N36-T3	CB,LL-KD	97	97	99	99	37	6,11*,2	* Pioneer 35A30		104	104	104	98	54	6*,12*,
NK Brand N38-B4	CB,LL,RR-KDH	97	97	99	37	6,11	6,11	* Pioneer 35A31	CB,LL-L	105	105	105	103	50	6*,12*,21*
NK Brand N39-K7	RR-H	98	98	96	37	11	11	* Pioneer 35D29	CB,LL-L	106	108		108	50	19*,21*
* NK Brand N44-M1	CB,LL,RR-KDH	101	101	101	100	37	6,12,23	* Pioneer 35F38		104	104	103	103	50	6,18*,21*
* NK Brand N45-A6	CB,LL-KD	101	101	101	37	6*,12	6*,12	Pioneer 35R58	CB-G	105	105		101	50	21
* NK Brand N45-T5		101	101		97	37	23*	* Pioneer 35Y61	CB,CR,LL,RR-LPJ	106	106	107	104	50	7*,12*,21
NK Brand N46-N9	CB,LL-KD	102	102	103	37	12	12	* Pioneer 35Y67	CB,LL-L	106	106		105	54	21*
NK Brand N47-F2	CB,LL-KD	102	102	103	37	12	12	* Pioneer 36H54	CB,LL-L	104	104	106		50	6,12*
NK Brand N48-R3	CB,LL-KD	103	103	102	37	6	6	* Pioneer 36W66		102	102		104	54	18*
* NK Brand N48-V8	CB,LL-KD	106	106		102	37	21*	* Pioneer 36W67	CB,LL-L	103	103	104		50	12*
* NK Brand N49-E3	Leafy-	97	106		103	37	18*,20*	* Pioneer 37A92	CB-G	97	97		92	65	24*
* NK Brand N51-T8	CB,LL-KD	104	104	103	104	37	6*,12*,21	Pioneer 37A93	CB,RR-GJ	97	97		92	50	24
* NK Brand N53-W3	RR-H	105	105	103	37	6*,12*	6*,12*	* Pioneer 37F74	CB,LL-L	100	100	99		50	11*
* NK Brand N58-D1	CB,LL-KD	107	107	104	37	12*,13*,14	12*,13*,14	* Pioneer 37Y12		97	97	95		50	16*
NK Brand N58-L8	RR-H	106	106	107	106	37	7,21	* Pioneer 38B85		96	96	95	95	54	11,16*,20*
NK Brand N65-C5	CB,LL-KD	109	109	109	37	7	7	* Pioneer 38B86	CB,LL,RR-LJ	98	98	96		50	11*,16*
NK Brand N65-Y3	CB,LL-KD	109	109		108	37	19,21	* Pioneer 38H65	CB,LL,RR-LJ	99	99	96	94	50	16*,23
* NK Brand N67-W5	CB,LL-KD	109	109	110	108	37	7*,19*	Pioneer 38K46		93	93	90		50	16,17
* NK Brand N72-G8		112		109	37	19*	19*	* Pioneer 38M59	CB-G	95	95	94	92	50	16,2
NK Brand N78-D6	Leafy,RR-H	116		112	37	19	19	* Pioneer 38W22	CB,LL-L	92	92	87	92	50	17*,
								Pioneer 39D82	CB-G	87	87	86		50	15,17
Nokomis AR21BL		108	110		3	10	10	Prairie Hybrids 0371		92	93		3	10	
North Gro 1541RR	RR-J	80		86	1	15	15	* Prairie Hybrids 1671		102	103		3	10*	
North Gro 1690RR	RR-J	84		86	1	15	15	* Prairie Hybrids 2431		103	101		3	10*	
* North Gro 5610RR	RR-J	100		101	98	1	11*,23	* Prairie Hybrids 3081		104	105		3	10*	
North Gro 5919RR	RR-J	104		106	1	21	21								
North Gro 5945YGPLRR	CB,CR,RR-GNJ	103	104		1	6,12	6,12	* Renk 6-354LLYGCB	CB,LL-K	88		92	48	24*	
								* Renk 6-662S		105		104	12	18*,21	

† Trait-Gene=Code: **bmr**-bm3=B; **IMI**:IT=C; **LL**-T25=D; **Bt-CB**-Bt176=F,Mon810=G,Bt11=K,TC1507=L; **RR**-MonGA21=H,Nk603=J, SytGA21=Q; **Leafy**-lfy=M; **Bt-CR**-Mon863=N, DAS591227=P; **Unknown**-X.

Table 2. Companies and hybrids included in the 2006 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes	Relative Maturity Seed				Tables	
		Co.	MN	GRM	SRM				Trt	Co.	MN	GRM		SRM
Renk RK244		83	83	85	12	17	Trelay 5N503	CB,RR-GJ	100	100	100	95	11,14	
* Renk RK282		88	85	89	48	22*	* Trelay 5N749	CB,RR-GJ	102	102	103	95	6*,12*	
* Renk RK288YGCB	CB-G	86	85	86	48	15,17*	Trelay 6K917	CB,CR,RR-GNJ	107	107	105	95	21	
Renk RK402		90	87		48	15	* Trelay 6N664	CB,RR-GJ	107	107	107	95	7,9*	
* Renk RK438RRYGCB	CB,RR-GJ	92	92	92	48	24*	Trelay 6R402	RR-J	106		104	95	21	
* Renk RK438YGCB	CB-G	92	92	91	48	16*,17*	* Trelay 7242	CB-G	90	90	88	95	15*	
* Renk RK452LLYGCB	CB,LL-K	94	94	92	48	24*	* Trelay 7560	CB-G	100	100	101	95	11,13*,14*	
* Renk RK488YGCB	CB-G	96	97	95	48	11,16*	Trelay 7N269	CB,RR-GJ	109	109		95	7	
* Renk RK575YGPL	CB,CR-GN	99	98	98	48	6*,11,20*	Trelay 7N412	CB,CR-GN	110	111		95	7	
* Renk RK632RRYGPL	CB,CR,RR-GNJ	102	100	102	101	48	6,12,21,	* Trelay 8K339	CB,CR-GN	113	110	111	95	7*,19
Renk RK636RRYGCB	CB,RR-GJ	103	103	104	48	6,12	Trelay 8R321	RR-J	114		110	95	19	
* Renk RK644YGCB	CB-G	105	103		48	6*,12*								
* Renk RK669		102	102	105	48	6,12*,18,21	UW EX27		112		111	1	19	
* Renk RK772RRYGPL	CB,CR,RR-GNJ	103	101		48	6*,12	* UW EX28		100		100	1	20*	
Renk RK789YGPL	CB,CR-GN	109	108	109	48	7	* UW EX30		106		110	1	19*	
Renk RK854		112	112	111	48	19	UW EX31		105		107	1	21	
* Renk RK888YGCB	CB-G	112		109	48	19*	UW EX32		100		100	1	20	
							* UW EX33		106		109	1	19*	
Spangler 150LG	CB,LL-K	82	85		37	15	UW EX34		101		105	1	21	
Spangler 1777		85	87		7	15	UW EX35		106		110	1	19	
* Spangler 324G	CB-G	92	94		37	11,16*								
* Spangler 407+	CB,CR-GN	99	96		52	11*	Unity Seeds 2200YGPL	CB,CR-GN	101	103		57	12	
* Spangler 457		103	104		14	12*	Unity Seeds 2210YGPL	CB,CR-GN	110	109		57	7	
* Spangler 458G	CB-G	104	103		37	6*,12*	Unity Seeds 6106C		106	108		95	7	
Spangler 667LG	CB,LL-K	107	110		37	7								
Spangler 686G	CB-G	109	110		94	7	Viking O.5305		103	103	99	3	10	
Spangler 7558G	CB-G	112	109		51	7	Viking O.7292		96	96		3	10	
Spangler RX217	RR-J	90	86		7	15								
* Trelay 2N553	CB,RR-GJ	90	90	90	95	22*	* Wolf River Valley 2490RRB	CB,RR-GJ	90	88		15	15*	
* Trelay 2R144	RR-J	90	90	91	95	22*	* Wolf River Valley 2797		97	95		15	16*	
* Trelay 4B268	CB-G	96	96	96	95	11,16*,								
* Trelay 4N627	CB,RR-GJ	98	98	99	95	11*	Wyffels W1959	CB,CR-GN	98	99		7	6	
* Trelay 5B353	CB-G	100	100	100	95	11*	Wyffels W2711	CB,CR,RR-GNJ	102	104		7	6,9	
* Trelay 5B423	CB-G	100	100	97	95	20*	Wyffels W3819	CB,CR-GN	104	100		7	6	
Trelay 5B739	CB-G	105	104	104	95	21	Wyffels W3945	CB,LL-L	105	105		7	6	
* Trelay 5K106	CB,CR,RR-GNJ	100	100	101	99	95	11,13*,20*	Wyffels W4167	CB,RR-GN	104	104		7	6,9

† Trait-Gene=Code: **bmr**-bm3=B; **IMI**:IT=C; **LL**-T25=D; **Bt-CB**-Bt176=F,Mon810=G,Bt11=K,TC1507=L; **RR**-MonGA21=H,Nk603=J, SytGA21=Q; **Leafy**-lfy=M; **Bt-CR**-Mon863=N, DAS591227=P; **Unknown**-X.