

Table 2. Companies and hybrids included in the 2008 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes†	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes†	Relative Maturity Seed				Tables
		Co.	GRM	SRM	Trt.				Co.	GRM	SRM	Trt.	
AgriGold A6225VT3	CB,CR,RR-GQJ	99	98		55	8	* Croplan Genetics 388TS	CB,CR,RR-GNJ	99	97		37	8*, 11*
* AgriGold A6279VT3	CB,CR,RR-GQJ	101	103		55	6, 9*	Croplan Genetics 4015VT3	CB,CR,RR-GQJ	98	96		37	8, 11
* AgriGold A6309BtRWR	CB,CR,RR-GNJ	103	103	101	55	6, 9, 15*	Croplan Genetics 421VT3	CB,CR,RR-GQJ	100	98		37	6, 8, 11
* AgriGold A6323CL	IMI-C	103	104	102	55	15*, 22*	Croplan Genetics 4801VT3	CB,CR,RR-GQJ	102	102		37	6
* AgriGold A6325VT3	CB,CR,RR-GQJ	104	104	101	55	6, 9, 15*	* Croplan Genetics 5338VT3	CB,CR,RR-GQJ	104	103	102	37	6*, 9, 15*, 18*
AgriGold A6394VT3	CB,CR,RR-GQJ	108	110	107	55	7, 13	* Croplan Genetics 5892VT3	CB,CR,RR-GQJ	107	108		37	7, 9, 25*
AgriGold A6399VT3	CB,CR,RR-GQJ	108	107		55	9	* Croplan Genetics 591VT3	CB,CR,RR-GQJ	108	107	107	37	7, 9, 16*, 18*
* AgriGold A6439VT3	CB,CR,RR-GQJ	109	109	107	55	7*, 13*, 16*	Croplan Genetics 6069AS3G	CB,CR,LL,RR-KRDS	110	109		37	14, 16, 18
AgriGold A6459VT3	CB,CR,RR-GQJ	110		112	55	14	* Croplan Genetics 6100VT2	CR,RR-QJ	110		108	37	14*, 16, 18*
AgriGold A6474VT3	CB,CR,RR-GQJ	111		111	55	14	Croplan Genetics 6226VT3	CB,CR,RR-GQJ	111	108		37	7
							* Croplan Genetics 645TS	CB,CR,RR-GNJ	111		110	37	14*, 16
* AgVenture AV7597V3R	CB,RR-GJ	106		106	118	13*	* Croplan Genetics 6818VT3	CB,CR,RR-GQJ	114	112	112	37	7, 14*
AgVenture AV7729CBLL	CB,LL-KD	107		107	118	13	Croplan Genetics 6831TS	CB,CR,RR-GNJ	112		112	37	14, 16
* AgVenture G6512	RR-J	97		103	118	13*, 15*	* Croplan Genetics 8221VT3	CB,CR,RR-GQJ	118		113	37	14*
* AgVenture R4926VBW	CB,CR,RR-GNJ	94		96	118	15*, 17*	* Croplan Genetics S4900	Leafy,CR,RR-MQJ	102		103	37	15*, 18*
* AgVenture R5677VBW	CB,CR,RR-GNJ	97		101	118	15*							
AgVenture R6608YBW	CB,RR-GJ	100		106	118	13	Crows 1466R	RR-J	85	86		37	10
* AgVenture R6864	RR-J	102		106	118	13*	* Crows 1685VT3	CB,CR,RR-GQJ	92	95		37	11*
* AgVenture RL5906HBW	CB,LL,RR-LDJ	97		99	118	15, 17*	* Crows 1725VT3	CB,CR,RR-GQJ	95	96		37	11*
							Crows 1807VT3	CB,CR,RR-GQJ	97	98		37	8
* Blue River 25M90		88	93		96	21*	* Crows 1928R	RR-J	99	98	96	37	8*, 11, 17*
Blue River 30B19		90	94		96	21	* Crows 2123VT3	CB,CR,RR-GQJ	101	102	104	37	6*, 18*
* Blue River 41R00		97		95	96	17*	* Crows 2155VT3	CB,CR,RR-GQJ	102	102		37	9*
* Blue River 42A32		96	96		96	21*	* Crows 3848VT3	CB,CR,RR-GQJ	105	103	105	37	6, 16*, 18*
Blue River 46L96		100		103	96	18	* Crows 4354VT3	CB,CR,RR-GQJ	107	109		37	7*
* Blue River 48B30		103	101	104	96	18*, 21*	Crows 4822B	CB-G	110	111	111	37	7, 14
* Blue River 56M30		106	104		96	21*	* Crows 4846T	CB,CR,RR-GNJ	110		109	37	16*
Blue River 59L36		108		107	96	18							
							Dahlco Seeds 6944VT3	CB,CR,RR-GQJ	95	95		55	8
Carharts Blue Top C9000Bt	CB-G	90	87		55	10, 12	Dahlco Seeds 7071VT3	CB,CR,RR-GQJ	107	109		55	7
* Carharts Blue Top CR1185V	CB,CR,RR-GQJ	85	87		2	10*, 12*	Dahlco Seeds 7810VT3	CB,CR,RR-GQJ	81	85		55	10
Carharts Blue Top CR1857R	CB,RR-GJ	90	95		55	8	* Dahlco Seeds 7871GTCBLL	CB,LL,RR-KDS	87	86		55	10*
Carharts Blue Top CR1995V	CB,CR,RR-GQJ	95	96		33	8, 11	* Dahlco Seeds 7881VT3	CB,CR,RR-GQJ	88	87		55	10*
Carharts Blue Top CR9680V	CB,CR,RR-GQJ	98	98		55	6, 8, 11	Dahlco Seeds 7931GTCBLL	CB,LL,RR-KDS	93	97		55	8
							* Dahlco Seeds 7951GTCBLL	CB,LL,RR-KDS	95	95		55	8*
* Cornelius C382YG	CB-G	103		105	99	13*	Dahlco Seeds 8001GTCBLL	CB,LL,RR-KDS	100	99		55	8
Cornelius C447VT3	CB,CR,RR-GQJ	104	103		99	6	Dahlco Seeds 8021VT3	CB,CR,RR-GQJ	102	102		55	9
* Cornelius C454XTLL	CB,CR,LL-LPD	105	106	107	99	6*, 13*	Dahlco Seeds 8811GTCBLL	CB,LL,RR-KDS	81	85		55	10
* Cornelius C462		104	104		99	22*	Dahlco Seeds 8841VT3	CB,CR,RR-GQJ	84	86		55	10
* Cornelius C466XTLL	CB,CR,LL-LPD	105	105		99	6*	Dahlco Seeds 8961VT3	CB,CR,RR-GQJ	96	97		55	8
Cornelius C547VT3	CB,CR,RR-GQJ	107	109		99	7							
* Cornelius C591		109	106		99	22*	* Dahlman D4356	CB,LL,RR-KDH	85	86	116		12*
* Cornelius C595YG	CB-G	109	109	105	99	7, 13*	Dahlman D4462	CB,LL,RR-KDH	88	86	116		12
Cornelius C649VT3	CB,CR,RR-GQJ	110	109	111	99	7, 14	Dahlman R4342CB	CB,RR-GJ	86	87	116		10
							Dahlman R4342VT3	CB,CR,RR-GNJ	86	87	99		10
Croplan Genetics 238VT3	CB,CR,RR-GQJ	85	85		37	10	Dahlman R4525VT3	CB,CR,RR-GNJ	90	89	99		10
Croplan Genetics 2520TS	CB,CR,RR-GNJ	86	86		37	10	Dahlman R45S15	CB,RR-GJ	90	89	116		10
Croplan Genetics 2538RR	RR-J	86	88		37	10	* Dahlman R4807VT3	CB,CR,RR-GNJ	95	95	99		11*
* Croplan Genetics 3114VT3	CB,CR,RR-GQJ	91	95	94	37	11*, 17*	* Dahlman R4928VT3	CB,CR,RR-GNJ	98	97	99		11*
* Croplan Genetics 3456VT3	CB,CR,RR-GQJ	95	97		37	11*							
* Croplan Genetics 3514VT3	CB,CR,RR-GQJ	95	98	95	37	11*, 17*	Dairyland HiDF3000-6	RR-J	100		103	33	15, 18
* Croplan Genetics 3724VT3	CB,CR,RR-GQJ	96	98		37	8, 11*	* Dairyland HiDF3007		107		110	33	16*

† Code = Trait(Gene); B = bmr(bm3); C = IMI(IT); D = LL(T25); F, G, K, L = Bt-CB: (Bt176, Mon810, Bt11, TC1507); H, J, S = RR(MonGA21, Nk603, SYTGA21); M = Leafy(lfy1); N, P, Q, R = Bt-CR (Mon863, DAS591227, Mon88017, MIR604); X = Unknown.

Table 2. Companies and hybrids included in the 2008 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes†	Relative Maturity			Seed Trt.	Tables	Brand Hybrid	Traits - Genes†	Relative Maturity			Seed Trt.	Tables
		Co. GRM	SRM	Trt.					Co. GRM	SRM	Trt.		
* Dairyland HiDF3008-4	CR,RR-NJ	108	106	33	13*	Dyna Gro CX08711	CB,CR,RR-GQJ	110	110	7	7		
Dairyland HiDF3012-6	RR-J	112	112	33	14								
* Dairyland HiDF3085-6	RR-J	85	86	33	19*	Farm Advantage 87A03	CB,CR,LL-KRD	103	103	43	6, 9		
* Dairyland HiDF3104		104	104	33	13*, 15*	Farm Advantage 87A10GL	CB,CR,LL,RR-KRDS	110	110	43	7		
* Dairyland HiDF3110-6	RR-J	110	112	33	14*	Farm Advantage 87A99GL	CB,CR,LL,RR-KRDS	99	99	43	6, 8		
Dairyland ST4006	CR,RR-NJ	106	107	33	7	* Farm Advantage 87X00	CB,CR,LL-LPD	100	102	43	6, 8*		
Dairyland ST7196	CB,RR-GJ	96	96	33	8, 11	* Farm Advantage 87X04	CB,CR,LL-LPD	104	106	43	6*		
* Dairyland ST7286	CB,RR-GJ	86	86	33	10, 12*								
* Dairyland ST7789	CB,RR-GJ	89	90	33	19*	Fielders Choice Direct NG	CB,CR,RR-GQJ	88	92	7	17		
Dairyland ST7891	CB,LL,RR-KDH	91	93	33	11	* Fielders Choice Direct NG	RR-J	107	106	7	13, 16*		
* Dairyland ST7985	CB,LL-KD	85	86	33	10, 12*	* Fielders Choice Direct NG	CB,CR,RR-GQJ	114	111	7	14, 16*		
* Dairyland ST8208	CB,CR,LL-LPD	108	110	33	7*	* Fielders Choice Direct NG	CB,CR,RR-GQJ	107	104	7	13*, 16, 18*		
Dairyland ST8404	CB,LL-LD	104	106	33	9								
* Dairyland ST9002	CB,LL-KD	102	101	33	6, 9*	Foundation Direct OR8511		104	106		21		
* Dairyland ST9005	CB,CR,RR-GNJ	105	103	33	6, 9*	Foundation Direct OR8511N		104	106		21		
Dairyland ST9006	CB,CR,RR-GQJ	106	104	33	9	* Foundation Direct OR8541N		105	105		21*		
* Dairyland ST9009	CB,CR,RR-GQJ	111	110	33	7*	* Foundation Direct OR8643		102	101		21*		
Dairyland ST9196	CB,CR,RR-GQJ	96	96	33	11	* Foundation Direct OR8643N		102	101		21*		
* Dairyland ST9789	CB,CR,RR-GQJ	89	88	33	10*, 12*	Foundation Direct OR8797		95	98		21		
* Dairyland ST9799	CB,CR,RR-GQJ	99	98	33	8, 11*, 25*	* Foundation Direct OR8800		95	95		21*		
* Dairyland ST9902	CB,RR-GJ	102	102	33	6*, 9*	* Foundation Direct OR8890		95	96		21*		
* Dekalb DKC33-54(RR2)	RR-J	83	85	7	12*	G2 Genetics 0A-016		116	111	7	14		
Dekalb DKC33-72(RR2)	RR-J	83	85	7	19	* G2 Genetics 1H-005AHxLL	CB,LL-LD	105	105	7	6*		
* Dekalb DKC35-19(RR2YGCB)	CB,RR-GJ	85	86	7	12*	G2 Genetics 1H-005HxLL	CB,LL-LD	105	105	7	6		
* Dekalb DKC38-89(VT3)	CB,CR,RR-GQJ	88	89	91	7, 10, 12*, 17*, 19*	G2 Genetics 1H-911HxLL	CB,LL-LD	111	109	7	7		
* Dekalb DKC41-60(VT3)	CB,CR,RR-GQJ	91	92	93	7, 11, 12*, 17*	G2 Genetics 1X-716HxTLL	CB,CR,LL-LPD	116	112	7	14		
* Dekalb DKC42-91(VT3)	CB,CR,RR-GQJ	92	92	93	7, 11*, 12, 17*, 19*	* G2 Genetics 1X-911HxTLL	CB,CR,LL-LPD	111	110	111	7, 14*		
Dekalb DKC43-27(VT3)	CB,CR,RR-GQJ	93	94	7	11	* G2 Genetics 3A-513RR	RR-J	113	110	7	7*		
* Dekalb DKC45-79(VT3)	CB,CR,RR-GQJ	95	97	98	7, 8, 11, 15*, 17	G2 Genetics 5H-212RRHx	CB,LL,RR-LDJ	112	111	7	7		
Dekalb DKC46-60(VT3)	CB,CR,RR-GQJ	96	96	7	8, 11	* G2 Genetics 5H-298RRHx	CB,LL,RR-LDJ	98	95	93	7, 11*, 17*		
* Dekalb DKC48-46(RR2YGPL)	CR,RR-GJ	98	100	99	7, 6, 13*, 17*	* G2 Genetics 5H-501RRHx	CB,LL,RR-LDJ	101	102	7	15*, 18*		
* Dekalb DKC50-44(VT3)	CB,CR,RR-GQJ	100	101	103	7, 6, 8, 13*, 15*, 20*, 23*, 24, 26*	G2 Genetics 5H-506RRHx	CB,LL,RR-LDJ	105	104	7	6		
						G2 Genetics 5H-508ARRHx	CB,LL,RR-LDJ	105	107	7	6		
* Dekalb DKC52-43(VT3)	CB,CR,RR-GQJ	102	100	7	6, 9, 20*	G2 Genetics 5H-508RRHx	CB,LL,RR-LDJ	108	110	7	7		
* Dekalb DKC52-59(VT3)	CB,CR,RR-GQJ	102	100	101	7, 6*, 9, 15*, 24*, 25*	* G2 Genetics 5H-702RRHx	CB,LL,RR-LDJ	102	103	7	6*		
* Dekalb DKC54-16(VT3)	CB,CR,RR-GQJ	104	104	103	7, 6, 9, 13, 15, 20*								
* Dekalb DKC57-79(RR2YGPL)	CB,CR,RR-GNJ	107	105	7	13*	Garst 82K79GT	RR-S	112	111	102	14		
* Dekalb DKC61-69(VT3)	CB,CR,RR-GQJ	111	109	111	7, 7*, 14*, 20*	* Garst 8384CBLLRW	CB,CR,LL-KRD	108	106	102	13*		
* Dekalb DKC63-42(VT3)	CB,CR,RR-GQJ	113	111	111	7, 7*, 14	Garst 8452CBLL	CB,LL-LD	107	108	102	13		
* Dekalb RX674VT3(VT3)	CB,CR,RR-GQJ	109	109	106	7, 7, 13*	* Garst 84H80-3000GT	CB,CR,LL,RR-KRDS	107	107	102	13*		
						Garst 85R29GTCBLL	CB,LL,RR-KDS	101	101	102	15		
* Dyna Gro 52P81	CB,RR-GJ	88	86	59	12*	* Garst 8687HXLLIT	CB,LL,IMI-LDC	100	102	102	15*		
* Dyna Gro 52V01	CB,CR,RR-GQJ	87	87	59	10, 12*	* Garst 87G94-3000GT	CB,CR,LL,RR-KRDS	94	95	102	17*		
* Dyna Gro 54T42	CB,RR-LJ	100	102	104	11, 18*	* Garst 88C97CBLL	CB,LL-LD	92	95	102	17*		
Dyna Gro 54V78	CB,CR,RR-GQJ	97	95	59	8, 11	* Garst 88F74GTRW	CR,RR-GJ	92	96	102	17*		
Dyna Gro 55V18	CB,CR,RR-GQJ	100	97	59	8, 11	* Garst 8907GTCBLL	CB,LL,RR-KDS	82	84	102	19*		
Dyna Gro 55V71	CB,CR,RR-GQJ	103	103	59	6, 9								
Dyna Gro 57P28	CB,RR-GJ	113	110	59	16	* Gold Country 10007VT3	CB,CR,RR-GQJ	100	97	7	8, 11*		
Dyna Gro 57V98	CB,CR,RR-GQJ	111	108	59	7	* Gold Country 10204VT3	CB,CR,RR-GQJ	102	101	7	9, 11*		
* Dyna Gro CX07102	CB,CR,RR-GNJ	102	103	7	18*	* Gold Country 9404VT3	CB,CR,RR-GQJ	94	97	116	8, 11*		
* Dyna Gro CX07305	CB,CR,LL-LPD	108	106	7	9*	* Gold Country 9608VT3	CB,CR,RR-GQJ	96	94	7	8*, 11*		

† Code = Trait(Gene); B = bmr(bm3); C = IMI(IT); D = LL(T25); F, G, K, L = Bt-CB: (Bt176, Mon810, Bt11, TC1507); H, J, S = RR(MonGA21, Nk603, SYTGA21); M = Leafy(lfy1); N, P, Q, R = Bt-CR (Mon863, DAS591227, Mon88017, MIR604); X = Unknown.

Table 2. Companies and hybrids included in the 2008 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes†	Relative Maturity			Seed Trt.	Tables	Brand Hybrid	Traits - Genes†	Relative Maturity			Seed Trt.	Tables
		Co.	GRM	SRM					Co.	GRM	SRM		
* Gold Country 9810VT3	CB,CR,RR-GQJ	98	97		7	8, 11*	Jung 7362VT3	CB,CR,RR-GQJ	94	97		55	8, 11
Great Lakes 3730G3VT3	CB,CR,RR-GQJ	87	87		20	10	* Jung 7426VT3	CB,CR,RR-GQJ	96	96		55	8*
* Great Lakes 4041G3VT3	CB,CR,RR-GQJ	90	89		20	10*	* Jung 7454VT3	CB,CR,RR-GQJ	98	98		55	8, 11*
Great Lakes 4415G3VT3	CB,CR,RR-GQJ	94	95		20	11	* Jung 7475VT3	CB,CR,RR-GQJ	100	100		55	11*
* Great Lakes 4689G3VT3	CB,CR,RR-GQJ	96	97	95	20	8, 17*	* Jung 7482VT3	CB,CR,RR-GQJ	100	97		55	11*
Great Lakes 4951BtRR	CB,RR-GJ	100	99		20	8	* Jung 7514VT3	CB,CR,RR-GQJ	102	101		55	11*
* Great Lakes 5335X2	CB,LL-LD	103	104	101	20	9*, 15*	Jung HDS6098QRWRR	CR,RR-QJ	98		99	55	15, 17
Great Lakes 5416G3VT3	CB,CR,RR-GQJ	104	103		20	6	Jung HDS65W44	CR,RR-QJ	94		96	55	17
Great Lakes 5995BtRW	CB,CR-GN	109	109		20	7	* Jung HDS7105VT3	CB,CR,RR-GQJ	104		104	55	13*, 15
							* Jung HDS76V78	CR,RR-QJ	107		107	55	13*
Growmark FS34S31	CB,RR-GJ	84	86		20	10	Kaltenberg K2403VT3	CB,CR,RR-GQJ	83	86		95	12
Growmark FS4189VT3	CB,CR,RR-GQJ	101	102		20	9	* Kaltenberg K2801RR	RR-J	84	85		95	12*
* Growmark FS44SV3	CB,CR,RR-GQJ	94	94	96	20	11*, 17*	* Kaltenberg K3023RRPlus	CB,CR,RR-GQJ	82	87		95	12*
* Growmark FS49SV3	CB,CR,RR-GQJ	99	97	95	20	11*, 17*	* Kaltenberg K3038LLBt11	CB,LL-KD	87	86		95	12*
* Growmark FS53SV3	CB,CR,RR-GQJ	103	103	101	20	9, 15*	Kaltenberg K3843RRPlus	CB,CR,RR-GQJ	95	96		95	8
* Growmark FS5559VT3	CB,CR,RR-GQJ	107	108	105	20	7, 13*	Kaltenberg K4013VT3	CB,CR,RR-GQJ	94	95		95	8, 11
Growmark FS55JV3	CB,CR,RR-GQJ	105	105		20	9	Kaltenberg K4263VT3	CB,CR,RR-GQJ	95	97		95	8, 11
* Growmark FS55SV3	CB,CR,RR-GQJ	105	104	106	20	9*, 16	Kaltenberg K4433VT3	CB,CR,RR-GQJ	100	97		95	8
* Growmark FS56SV3	CB,CR,RR-GQJ	106	105	107	20	9*, 16	Kaltenberg K5163VT3	CB,CR,RR-GQJ	103	103		95	9
* Growmark FS58SV3	CB,CR,RR-GQJ	108	109	106	20	7, 13*	Kaltenberg K5232RRLlBtHx	CB,LL,RR-LDJ	103	104		95	9
* Growmark FS60AV3	CB,CR,RR-GQJ	110	109	111	20	7, 14*	Kaltenberg K5578LLBt11	CB,LL-KD	105	106		95	9
* Growmark FS60JV3	CB,CR,RR-GQJ	110	109	111	20	7*, 14	Kaltenberg K5823VT3	CB,CR,RR-GQJ	106	108		95	7
Growmark FS61AV3	CB,CR,RR-GQJ	111		112	20	14	Kaltenberg K6165RRLlHXT	CB,CR,LL,RR-LPDJ	108	109		95	7
Growmark FS6299VT3	CB,CR,RR-GQJ	110	109	111	20	7, 14	* Kaltenberg K6355RRLlBtHx	CB,LL,RR-LDJ	109	110	107	95	7, 16*
* Growmark FSE3802	CB,LL,RR-LDJ	96	95	94	20	11, 17*	Kaltenberg K6648LLBt11	CB,LL-KD	110	110		95	7
							Kaltenberg K8095LFRR	Leafy,RR-MJ	95		95	95	19
* Hughes 2823HXRR	CB,LL,RR-LDJ	101	102		37	6*	* Kaltenberg K8099LFRR	Leafy,RR-MJ	99		95	95	19*
Hughes 3208CBLL	CB,LL-LD	102	102		37	6	Kaltenberg K8105LFRR	Leafy,RR-MJ	105		106	95	18
Hughes 3878CBLLRW	CB,CR,LL-KRD	104	104		37	6	Kaltenberg K8112LFRR	Leafy,RR-MJ	112		109	95	16
* Hughes 4592VT3	CB,CR,RR-GQJ	105	105	105	37	6, 13*							
* Hughes 4898HXT	CB,CR,LL-LPDJ	106	109		37	7*	Kruger K1086RR	RR-J	86	84		37	12
Hughes 6347VT3	CB,CR,RR-GQJ	110	110		37	7	Kruger K1087RR	RR-J	87	86		37	10, 12
							Kruger K1093RR	RR-J	93	89		37	12
ISU AR16026:S1704OP		108	109			21	Kruger K1178RR	RR-J	78	85		37	12
Johnson Seeds 2195VT3	CB,CR,RR-GQJ	90	89		55	10, 12	* Kruger K1295RR	RR-J	95	96		37	8, 11*, 26
Johnson Seeds 2241RR	RR-J	97	95		55	11	Kruger K1490RR	RR-J	90	87		37	10, 12
Johnson Seeds 7100CBLLRW	CB,CR,LL-KRD	100	98		55	8, 11	* Kruger K1500RR	RR-J	99	97		37	6*, 8*, 11*, 23*, 26
Johnson Seeds 7189CBLLRW	CB,CR,LL-KRD	89	88		55	10, 12	Kruger K1584RR	RR-J	84	86		37	10, 12
* Johnson Seeds 7192GTCBLL	CB,LL,RR-KDS	92	91		55	11, 12*	Kruger K1780RR	RR-J	80	85		37	12
Johnson Seeds 7199CBLLRW	CB,CR,LL-KRD	99	98		55	8, 11	* Kruger K2086RRYGCB	CB,RR-GJ	86	85		37	10, 12
* Johnson Seeds 7804RRCRW	CR,RR-NJ	104	100		55	9, 11, 25*	* Kruger K2087RRYGCB	CB,RR-GJ	87	86		37	10, 12*
Johnson Seeds 7881GTCBLL	CB,LL,RR-KDS	81	85		55	12	* Kruger K2090RRYGCB	CB,RR-GJ	90	86		37	10*, 12
* Johnson Seeds 7886VT3	CB,CR,RR-GQJ	86	86		55	10, 12*	* Kruger K3300RRHx	CB,RR-LJ	100	102		37	6*, 8*
Johnson Seeds 7892CBLLGT	CB,LL,RR-KDS	93	94		55	11	* Kruger K5388YGCB	CB,RR-GJ	88	85		37	10, 12*
Johnson Seeds 7896VT3	CB,CR,RR-GQJ	96	97		55	11	* Kruger K6006VT3	CB,CR,RR-GQJ	106	104		37	7, 9*, 20*, 23*, 24, 26
							* Kruger K6007TS	CB,CR,RR-GNJ	107	107		37	7, 20*, 23*
Jung 3385RR	RR-J	94	95		55	8, 11	* Kruger K6008VT3	CB,CR,RR-GQJ	108	107		37	7, 20*, 23*
Jung 4183RRYGCB	CB,RR-GJ	84	86		55	10	* Kruger K6011TS	CB,CR,RR-GNJ	111	108		37	7, 20, 23*
Jung 7171VT3	CB,CR,RR-GQJ	83	87		55	10	* Kruger K6093VT3	CB,CR,RR-GQJ	93	97		37	8, 11, 24*, 26
Jung 7344VT3	CB,CR,RR-GQJ	91	96		55	11	* Kruger K6094VT3	CB,CR,RR-GQJ	94	97		37	8, 11*, 24, 26

† Code = Trait(Gene); B = bmr(bm3); C = IMI(IT); D = LL(T25); F, G, K, L = Bt-CB: (Bt176, Mon810, Bt11, TC1507); H, J, S = RR(MonGA21, Nk603, SYTGA21); M = Leafy(lfy1); N, P, Q, R = Bt-CR (Mon863, DAS591227, Mon88017, MIR604); X = Unknown.

Table 2. Companies and hybrids included in the 2008 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes†	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes†	Relative Maturity Seed				Tables
		Co.	GRM	SRM	Trt.				Co.	GRM	SRM	Trt.	
* Kruger K6097VT3	CB,CR,RR-GQJ	97	97	37	8*, 11*, 24*, 26*		Lemke 2183RRBt	CB,RR-GJ	87	87	95	10	
* Kruger K6102VT3	CB,CR,RR-GQJ	102	102	37	6, 9, 20*, 23*, 24, 26*		* Lemke 3087VT3	CB,CR,RR-GQJ	90	89	95	10*	
* Kruger K6107VT3	CB,CR,RR-GQJ	107	106	37	7, 20*, 23*		* Lemke 4097VT3	CB,CR,LL,RR-GQDJ	95	95	95	11*	
* Kruger K6111TS	CB,CR,RR-GNJ	111	108	37	7, 20*, 23*		* Lemke 5141Hx	CB,LL-LD	98	99	95	8*	
* Kruger K6208VT3	CB,CR,RR-GQJ	108	108	37	7*, 20*, 23*		Lemke 5160		100	99	95	8	
* Kruger K6210TS	CB,CR,RR-GNJ	110	110	37	7*, 20*, 23*		* Lemke 6110		105		106	95	
* Kruger K6212TS	CB,CR,RR-GNJ	112	112	37	7, 20*, 23*		Lemke 6187VT3	CB,CR,RR-GQJ	105	105	95	9	
Kruger K6298VT3	CB,CR,RR-GQJ	98	98	37	6, 8, 11, 24, 26		Lemke 7068Bt	CB-G	112		109	95	
* Kruger K6400TS	CB,CR,RR-GNJ	100	100	37	6, 8, 11, 20*, 23*, 24, 26		Lemke 7158-3000	CB,CR,LL,RR-KRDS	110	106	95	9	
* Kruger K6401VT3	CB,CR,RR-GQJ	101	103	37	6, 9*, 20*, 23*, 24*, 26*		LG Seeds LG2358VT3	CB,CR,RR-GQJ	85	86	54	10	
* Kruger K6411VT3	CB,CR,RR-GQJ	111	111	37	7, 20*, 23*		LG Seeds LG2362VT3	CB,CR,RR-GQJ	85	86	54	10	
* Kruger K6412VT3	CB,CR,RR-GQJ	112	110	37	7, 20*, 23*		* LG Seeds LG2411VT3	CB,CR,RR-GQJ	91	95	54	11*	
* Kruger K6499VT3	CB,CR,RR-GQJ	99	99	37	6, 8, 11*, 20*, 23*, 24*, 25*, 26		* LG Seeds LG2426VT3	CB,CR,RR-GQJ	94	95	54	11*	
Kruger K6503TS	CB,CR,RR-GNJ	103	102	37	6, 9		LG Seeds LG2469VT3	CB,CR,RR-GQJ	97	97	54	8	
* Kruger K6606VT3	CB,CR,RR-GQJ	106	105	37	7, 9*, 20*, 23*, 24, 25*, 26		LG Seeds LG2477VT3	CB,CR,RR-GQJ	99	99	54	8	
Kruger K6697VT3	CB,CR,RR-GQJ	97	96	37	8, 11, 24, 26		* LG Seeds LG2496VT3	CB,CR,RR-GQJ	99	98	54	8*	
* Kruger K7010YG+	CB,CR-GN	110	109	37	7, 20*		* LG Seeds LG2509CL	IMI-C	104		101	54	
Kussmaul GL790GT	RR-H	90	92	64	12		* LG Seeds LG2514VT3	CB,CR,RR-GQJ	104	104	54	9*	
* Kussmaul GL890Triple	CB,CR,LL,RR-KNDJ	90	87	64	12*		LG Seeds LG2532VT3	CB,CR,RR-GQJ	105	103	59	6	
Kussmaul GL895Triple	CB,CR,RR-KRDS	95	95	64	11		* LG Seeds LG2545VT3	CB,CR,RR-GQJ	107	110	54	7*	
* Kussmaul GL899GT	RR-J	99	99	64	8*		LG Seeds LG2548		108	111	54	7	
Kussmaul K700		100	99	68	8		LG Seeds LG2552VT3	CB,CR,RR-GQJ	110		112	54	
Kussmaul K710		110	108	68	7				105	106	7	22	
Kussmaul K795		95	95	68	11		* Masters Choice MC530		92	87	7	19*	
Kussmaul K-HC2		108	109	68	7		* Masters Choice SP490VT3	CB,CR,RR-GQJ	110		111	7	
* Kussmaul KX-706Hx	CB,LL-LD	105	106	64	6*		Masters Choice SP573VT3	CB,CR,RR-GQJ	110		111	7	
Kussmaul RFS610RR	RR-J	110		112	68				105	103		21	
Kussmaul SB590VTTriple	CB,CR,RR-GQJ	90	90	64	10		Michael Fields GEM2xCH5N2		105	102		21	
Kussmaul SB605VT3	CB,CR,RR-GQJ	106	107	64	7		Michael Fields M632xM73xC		105	104		21	
Kussmaul SB710RR	RR-J	110		112	68		Michael Fields M73-2xCH5N		105	104		21	
Legacy Seeds L2820CBLLRW	CB,CR,LL-KRD	86	88	37	12		Midwest Genetics 68101R	RR-J	85	86	37	12	
* Legacy Seeds L2850VT3	CB,CR,RR-GQJ	88	87	7	12*		* Midwest Genetics 69205VT3	CB,CR,RR-GQJ	92	92	37	11, 12*	
* Legacy Seeds L2927VT3	CB,CR,RR-GQJ	92	95	91	7	11, 19*	* Midwest Genetics 69575VT3	CB,CR,RR-GQJ	95		91	37	
* Legacy Seeds L3295VT3	CB,CR,RR-GQJ	96	97	7	6, 8, 11*		* Midwest Genetics 69704VT3	CB,CR,RR-GQJ	97	98	96	37	
* Legacy Seeds L3333VT3	CB,CR,RR-GQJ	95	95	7	11*		* Midwest Genetics 70006R	RR-J	99	98	97	37	
Legacy Seeds L3750VT3	CB,CR,RR-GQJ	97	98	7	8, 11		* Midwest Genetics 70505VT3	CB,CR,RR-GQJ	101	104	37	8*, 11, 19, 25*	
Legacy Seeds L4050VT3	CB,CR,RR-GQJ	100	98	7	8		* Midwest Genetics 72116VT3	CB,CR,RR-GQJ	102	102	102	37	
Legacy Seeds L4545LFY	Leafy-M	100		103	18		* Midwest Genetics 76126VT3	CB,CR,RR-GQJ	105	103	106	37	
* Legacy Seeds L5350CBLLGT	CB,LL,RR-KDH	104	104	103	37	6*, 9*, 15*, 18*	* Midwest Genetics 76485VT3	CB,CR,RR-GQJ	107	109	105	37	
Legacy Seeds L5950VT3	CB,CR,RR-GQJ	107	106	7	7, 9		Midwest Genetics 76804Y	CB,CR-GN	109	108	108	37	
* Legacy Seeds L6600Hx	CB,LL-LD	110		110	37	16*	* Midwest Genetics 76996VT3	CB,CR,RR-GQJ	109	108	105	37	
* Legend Seeds LR9707RRHxT	CR,LL,RR-PDJ	107	105	104	95	11, 18*	* Midwest Genetics 77125T	CB,CR,RR-GNJ	110	111	37	7, 13*	
* Legend Seeds LR9796VT3	CB,CR,RR-GQJ	96	97	94	95	11, 17*			100	99	37	8, 11	
Legend Seeds LR9798VT3	CB,CR,RR-GQJ	98	97	96	95	11, 17*	* Mycogen 2A517	CB,LL-LD	107	108	37	7*	
* Legend Seeds LS9703VT3	CB,CR,RR-GQJ	103	102	103	95	11, 18*	* Mycogen 2C598	CB,CR,LL,RR-LPDH	92	95	37	11	
							Mycogen 2D326	CB,CR,RR-GNJ	102	101	37	6, 9	
							Mycogen 2D519	CB,CR,RR-GQJ	92	91	37	11*, 12*	
							* Mycogen 2J337	CB,CR,RR-GQJ	96	95	37	8, 11	
							Mycogen 2J463	CB,CR,RR-GQJ					

† Code = Trait(Gene); B = bmr(bm3); C = IMI(IT); D = LL(T25); F, G, K, L = Bt-CB: (Bt176, Mon810, Bt11, TC1507); H, J, S = RR(MonGA21, Nk603, SYTGA21); M = Leafy(lfy1); N, P, Q, R = Bt-CR (Mon863, DAS591227, Mon88017, MIR604); X = Unknown.

Table 2. Companies and hybrids included in the 2008 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes†	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes†	Relative Maturity Seed				Tables
		Co.	GRM	SRM	Trt.				Co.	GRM	SRM	Trt.	
Mycogen 2K152	RR-H	82	87		37	10, 12	* NuTech 3T-096VT3	CB,CR,RR-GQJ	96	97	96	102	11, 17*
* Mycogen 2K660	CB-G	107	109		37	7*	* NuTech 3T-098AVT3	CB,CR,RR-GQJ	98	97		102	11*
* Mycogen 2M495	CB,CR,RR-GQJ	99	97		37	8, 11*	* NuTech 3T-098VT3	CB,CR,RR-GQJ	98	98	95	102	11, 17*
* Mycogen 2P172	RR-J	84	86		37	10, 12*	NuTech 3T-109VT3	CB,CR,RR-GQJ	109	109		7	7
* Mycogen 2P483	CB,CR,RR-GNJ	97	95		37	8, 11*	* NuTech 3T-110VT3	CB,CR,RR-GQJ	110	109		7	7*
Mycogen 2P535	CB,CR,LL,RR-LPDJ	103	103		37	6, 9	NuTech 3T-310VT3	CB,CR,RR-GQJ	110	109	112	7	7, 14
Mycogen 2R428	CB,CR,RR-GNJ	96	97		37	8	NuTech 3T-311VT3	CB,CR,RR-GQJ	111	109	111	7	7, 14
Mycogen 2R693	CB,CR,LL,RR-LPDJ	109	110		37	7	NuTech 3T-388VT3	CB,CR,RR-GQJ	88	86		7	10
* Mycogen F2F449	BMR,CB,LL,RR-BLDH	104		107	37	18*	* NuTech 3T-393VT3	CB,CR,RR-GQJ	90	89		7	10*, 12
* Mycogen F2F487	BMR,CB,LL-BLD	99		99	37	17*	NuTech 3T-415VT3	CB,CR,RR-GQJ	115	112		7	7
* Mycogen F2F633	BMR,CB,LL-BLD	110		110	37	16*, 18	NuTech 3T-484VT3	CB,CR,RR-GQJ	84	87		7	12
* Mycogen F2F699	BMR,CB,LL-BLD	113		113	37	16*	NuTech 3T-500VT3	CB,CR,RR-GQJ	100	99		7	11
* Mycogen F2F725	BMR,CB,LL-BLD	113		113	37	14*	NuTech 3T-510VT3	CB,CR,RR-GQJ	110	108		7	7
Mycogen F2F797	BMR-B	115		114	37	14	* NuTech 3T-514VT3	CB,CR,RR-GQJ	114	111	112	7	7*, 14
* Mycogen T29516	CB,CR,RR-GNJ	102		101	37	15*	NuTech 3T-710VT3	CB,CR,RR-GQJ	110	111		7	7
* Mycogen TMF2L412		92		94	37	17*, 19*	NuTech 3T-808VT3	CB,CR,RR-GQJ	108	109		102	7
Mycogen TMF2N422	RR-H	95		94	37	17	NuTech 3T-912VT3	CB,CR,RR-GQJ	112	110		7	7
Mycogen TMF2P719	CB,CR,LL-LPD	110		112	37	14	NuTech 3W-403RRYGRW	CR,RR-NJ	103	102		7	6
* Mycogen TMF2Q716	CB,CR,LL,RR-LPDJ	110		109	37	14*, 16	NuTech 5N-898GTCBLLRW	CB,CR,LL,RR-KRDS	98	98		7	11
* Mycogen TMF2Q733	CB,CR,LL-LPD	113		112	37	14*	NuTech 5X-008RRHxT	CB,CR,LL,RR-LPDJ	108	110		7	7
Mycogen TMF2Q759	CB,CR,LL,RR-LPDH	113		113	37	14	* NuTech 5X-201+HxTRR	CB,CR,LL,RR-LPDJ	101		101	7	15*
* Mycogen TMF2W726	CB,CR,LL-LPD	111		113	37	14*	NuTech 5X-512RRHxT	CB,CR,LL,RR-LPDJ	112	111	111	7	7, 14
* Mycogen X28566	CB,CR,LL-LPD	105	105	105	37	6, 9, 16*, 18	OBrien OB1028HxTLL	CB,CR,LL-LPD	100	101		54	8
* Mycogen X29209	CB,CR,LL-LPD	86	87		37	10, 12*	OBrien OB1059CBLLRWGT	CB,CR,LL,RR-KRDS	105	105		117	6
* Next Generation Seed 3787	CB,LL,RR-KDH	87	88	86	7	12*	OBrien OB1060HxLL	CB,LL-LD	105	104		64	6
Next Generation Seed 4899		99	100		37	8, 22	OBrien OB1087		108	109		43	22
Next Generation Seed 5802	CB,LL,RR-KDR	102	103		7	6	OBrien OB966		96	98		55	8
* Next Generation Seed 5804	CB,LL,RR-KDH	104		103	37	18*	OBrien OBX1009Hx	CB,LL-LD	100	99		117	8
Next Generation Seed 5806	CB,CR,LL,RR-KRDS	106	109		106	7, 20	* Organic UTC		103	98			21*
* Next Generation Seed 5808		108		106	7	13*, 16*	* Organic UTC-Hand Weed		103	99			21*
* NuTech 0C-404AYGCB	CB-G	104	103		7	6*	Pioneer 32T84	CB,LL,RR-LDJ	115		112	7	14
* NuTech 0C-404YGCB	CB-G	104	103		7	6*	* Pioneer 34A89	CB,CR,LL,RR-LPDJ	109		109	7	16*
NuTech 0C-413YGCB	CB-G	113	110		7	7	* Pioneer 35F40	CB,LL,RR-LDJ	104	105	103	7	6*, 9*, 15*, 18*, 20*, 24, 25*
* NuTech 1B-887CBLL	CB,LL-KD	85	86		7	12*	* Pioneer 36V53	CB,LL,RR-LDJ	102	103		7	9*
NuTech 1N-398CBLLRW	CB,CR,LL-KRD	98	98		102	11	* Pioneer 37Y14	CB,CR,LL,RR-LPDJ	99	97		7	8, 11*
NuTech 1N-887CBLLRW	CB,CR,LL-KRD	87	88		7	12	* Pioneer 38N87	CB,LL,RR-LDJ	92	91		7	11*, 12
* NuTech 3A-093RR	RR-J	90	89		7	10, 12*	Pioneer 39D82	CB-G	87	87		7	12
NuTech 3A-095RR	RR-J	95	96		7	11							
NuTech 3A-313RR	RR-J	113		112	102	14	* Prairie Hybrids 0336		96	95			21*
* NuTech 3A-390RR	RR-J	90	91		7	10, 12, 17*, 19*	Prairie Hybrids 2241		102	101			21
* NuTech 3A-804GT	RR-H	104		104	102	13*	* Prairie Hybrids 3081		100	104			21*
NuTech 3C-104RRYGCB	CB,RR-GJ	104	105		7	6	Prairie Hybrids 541		93	95			21
* NuTech 3C-300RRYGCB	CB,RR-GJ	100	99		7	11*							
* NuTech 3C-408RRYGCB	CB,RR-GJ	108	109		7	7*							
* NuTech 3P-191RRYGPL	CB,CR,RR-GNJ	90	90	90	102	10*, 12, 19*	* Renk RK268RRYGPL	CB,CR,RR-GNJ	84	87	88	55	12*, 19*
* NuTech 3P-302RRYGPL	CB,CR,RR-GQJ	102		101	102	15*	* Renk RK292CBLL	CB,LL-KD	85	87	86	55	10, 12, 19*
* NuTech 3P-616RRYGPL	CB,CR,RR-GNJ	116		112	102	14*	Renk RK501YGCB	CB-G	95	97	93	55	8, 11, 17
NuTech 3P-708BRRYGPL	CB,CR,RR-GNJ	105	106		7	6	* Renk RK570VT3	CB,CR,RR-GQJ	95	95		55	11*
NuTech 3P-708RRYGPL	CB,CR,RR-GNJ	108	110		7	7	Renk RK584CBLL	CB,LL-KD	99	99		55	6, 8, 11

† Code = Trait(Gene); B = bmr(bm3); C = IMI(IT); D = LL(T25); F, G, K, L = Bt-CB: (Bt176, Mon810, Bt11, TC1507); H, J, S = RR(MonGA21, Nk603, SYTGA21); M = Leafy(lfy1); N, P, Q, R = Bt-CR (Mon863, DAS591227, Mon88017, MIR604); X = Unknown.

Table 2. Companies and hybrids included in the 2008 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

Brand Hybrid	Traits - Genes†	Relative Maturity Seed				Tables	Brand Hybrid	Traits - Genes†	Relative Maturity Seed				Tables
		Co.	GRM	SRM	Trt.				Co.	GRM	SRM	Trt.	
* Renk RK616VT3	CB,CR,RR-GQJ	101	100	103	55	6, 9, 11, 18*	Viking O.5305		103	103			21
Renk RK670VT3	CB,CR,RR-GQJ	103	103		55	6, 9	* Viking O.6710		98	96			21*
Renk RK686VT3	CB,CR,RR-GQJ	102	104		55	6, 9	Viking O.7120		95	98			21
* Renk RK692CBLLRW	CB,CR,LL-KRD	105	105	105	55	6, 9, 16*, 18*, 25*	* Viking O.7292		96	95			21*
* Renk RK698RRYGRW	CR,RR-NJ	104	102	102	55	6, 9, 15*, 18*							
Renk RK719VT3	CB,CR,RR-GQJ	107	108		55	7	Wolf River Valley 2982	CB,RR-GJ	82	88		15	10
* Renk RK760RRYGCB	CB,RR-GJ	106	107	108	55	7, 16*	Wolf River Valley 2987	CB,LL,RR-KDH	87	86		15	10
* Renk RK770VT3	CB,CR,RR-GQJ	108	110	108	55	7, 13, 16*	Wolf River Valley 2995	CB,LL,RR-KDH	95	95		15	11
Renk RK829VT3	CB,CR,RR-GQJ	112	110	112	55	7, 14							
* Renk RK844VT3	CB,CR,RR-GQJ	112	110	112	55	7*, 14	Wyffels W1941	CB,CR,RR-GQJ	98	101		7	6
							Wyffels W2329	CB,CR,LL,RR-KRDH	100	100		102	6
* Spangler LFT49	Leafy-M	102		103	118	13*	Wyffels W2681	CB,CR,RR-GQJ	101	102		7	6
							Wyffels W5159	CB,CR,RR-KRDH	108	109		102	7
* Trelay 4T105	CB,CR,RR-GQJ	97	97		95	11*	Wyffels W5281	CB,CR,RR-GQJ	108	107		7	7
* Trelay 4T722	CB,CR,RR-GQJ	97	99		95	8, 11*, 25*	Wyffels W5641	CB,CR,RR-GQJ	109	109		7	7
* Trelay 5T128	CB,CR,RR-GQJ	100	98		95	8*	Wyffels W7251	CB,CR,RR-GQJ	112	111		7	7
* Trelay 5T429	CB,CR,RR-GQJ	101	103		95	9*							
Trelay 5T750	CB,CR,RR-GQJ	99	97		95	8	YIELDirect 3L73		101	101		54	22
* Trelay 6T226	CB,CR,RR-GQJ	106		105	95	16*	* YIELDirect 4L15		105	105		54	22*
Trelay 6T672	CB,CR,RR-GQJ	107		107	95	16	YIELDirect 4L92		107	108		54	22
Trelay 7K456	CB,CR,RR-GNJ	110		111	95	14	YIELDirect 4M62		104	106		54	22
Trelay 7T231	CB,CR,RR-GQJ	111		112	95	14	YIELDirect 5L96		109	108		54	22
Trelay 8T339	CB,CR,RR-GQJ	113	108		95	7	* YIELDirect 6E11		109	106		54	22*
UW EX27		112		112	1	14							
* UW EX31		110		111	1	16*							
UW EX35		110		112	1	14							
* Viking O.4700		106	103			21*							

Table 3. List of seed treatments used on corn hybrids entered in the 2008 UW corn trials.

Seed Trt.	CommonName	Seed Trt.	CommonName
1	Captan	64	ApronXL+Lorsban+Maxim+Poncho250
2	Maxim	68	ApronXL+Lorsban+Maxim
7	Poncho250	95	Captan+Poncho250
15	Actellic+MaximXL	96	Natural II
20	ApronFL+Maxim+Poncho250+Trilex	99	MaximXL+Poncho250+Trilex
33	Allegiance+Captan+Poncho250	102	Cruiser250
37	ApronXL+Cruiser+Dynasty+MaximXL	106	Actellic+Allegiance+Captan+Trilex
43	Cruiser+MaximXL	116	MaximXL+Trilex
54	ApronXL+Maxim	117	Actellic+MaximXL+Poncho250+Trilex
55	ApronXL+Maxim+Poncho250	118	ApronXL+Dynasty+Lorsban+Maxim+Poncho250+Trilex
59	Poncho250+Trilex		

† Code = Trait(Gene); B = bmr(bm3); C = IMI(IT); D = LL(T25); F, G, K, L = Bt-CB: (Bt176, Mon810, Bt11, TC1507); H, J, S = RR(MonGA21, Nk603, SYTGA21); M = Leafy(lfy1); N, P, Q, R = Bt-CR (Mon863, DAS591227, Mon88017, MIR604); X = Unknown.