

Table 2. Companies and hybrids included in the 2009 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

| Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables | Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables |
|----------------------------|------------------|------------------------|-----|-----|------|-------------|----------------------------|------------------|------------------------|-----|-----|------|-------------|
| | | Co. | GRM | SRM | Trt. | | | | Co. | GRM | SRM | Trt. | |
| AgPro 33C08 | | 108 | 106 | | 54 | 20 | | | | | | | |
| *AgVenture G6512 | RR-J | 100 | | 104 | 99 | 13*,15 | *Cornelius C286-300GT | CB,CR,LL,RR-KRDH | 100 | 98 | | 114 | 8* |
| *AgVenture GL7729ABW | CB,CR,LL,RR-KRDS | 107 | | 107 | 99 | 13* | *Cornelius C317XTLLRR | CB,CR,LL,RR-LPDJ | 100 | 100 | | 99 | 8* |
| *AgVenture R3530VBW | CB,CR,RR-GQJ | 88 | | 96 | 99 | 15,17* | *Cornelius C319VT3 | CB,CR,RR-GQJ | 101 | 106 | | 99 | 9 |
| *AgVenture R6864 | RR-J | 103 | | 105 | 99 | 13*,15* | *Cornelius C428XTLLRR | CB,CR,LL,RR-LPDJ | 104 | 104 | 106 | 99 | 6,9,13 |
| AgVenture RL5381HB | CB,LL,RR-LDJ | 95 | | 101 | 37 | 15 | *Cornelius C447VT3 | CB,CR,RR-GQJ | 105 | 104 | | 99 | 6 |
| AgriGold A6220VT3 | CB,CR,RR-GQJ | 98 | 99 | | 99 | 6,8 | *Cornelius C454XTLL | CB,CR,LL,RR-LPDJ | 105 | 104 | 107 | 99 | 6,13* |
| AgriGold A6225VT3 | CB,CR,RR-GQJ | 98 | 100 | | 99 | 6,8 | *Cornelius C462 | | 104 | 105 | | 99 | 20* |
| AgriGold A6279VT3 | CB,CR,RR-GQJ | 100 | 101 | | 99 | 6,8 | *Cornelius C547VT3 | CB,CR,RR-GQJ | 107 | 109 | | 99 | 7* |
| *AgriGold A6309VT3 | CB,CR,RR-GQJ | 103 | 103 | | 99 | 6*,9* | *Cornelius C591 | | 109 | 106 | | 99 | 20* |
| AgriGold A6320VT3 | CB,CR,RR-GQJ | 103 | 103 | | 99 | 6,9 | *Cornelius C595YG | CB-G | 108 | 109 | 106 | 99 | 7,13* |
| *AgriGold A6323CL | C | 103 | | 105 | 99 | 13*,15 | *Cornelius C649VT3 | CB,CR,RR-GQJ | 111 | 109 | | 99 | 7* |
| AgriGold A6325VT3 | CB,CR,RR-GQJ | 105 | 104 | | 99 | 6,9 | *Croplan Genetics 238VT3 | CB,CR,RR-GQJ | 85 | 85 | | 102 | 12* |
| AgriGold A6399VT3 | CB,CR,RR-GQJ | 108 | 108 | | 99 | 7,9 | *Croplan Genetics 2520VT3 | CB,CR,RR-GQJ | 86 | 87 | | 102 | 12 |
| *AgriGold A6439VT3 | CB,CR,RR-GQJ | 109 | | 107 | 99 | 13*,16 | *Croplan Genetics 3114VT3 | CB,CR,RR-GQJ | 91 | 90 | 95 | 102 | 11,12*,17 |
| *AgriGold A6459VT3 | CB,CR,RR-GQJ | 110 | | 111 | 99 | 14*,16* | *Croplan Genetics 3456VT3 | CB,CR,RR-GQJ | 95 | 97 | | 102 | 11 |
| *AgriGold A6533VT3 | CB,CR,RR-GQJ | 113 | | 111 | 99 | 14*,16 | *Croplan Genetics 3514VT3 | CB,CR,RR-GQJ | 95 | 96 | 95 | 102 | 11,17* |
| *Blue River 19K19 | | 84 | 85 | | | 23* | *Croplan Genetics 3724VT3 | CB,CR,RR-GQJ | 96 | 97 | | 102 | 11 |
| *Blue River 25M90 | | 86 | 90 | | | 23* | *Croplan Genetics 388TS | CB,CR,RR-GNJ | 99 | 97 | | 102 | 8,11 |
| *Blue River 30A12 | | 90 | 85 | | | 23* | *Croplan Genetics 421VT3 | CB,CR,RR-GQJ | 100 | 98 | | 102 | 8,11 |
| Blue River 33L90 | | 93 | | 91 | 96 | 19 | *Croplan Genetics 4338VT3 | CB,CR,RR-GQJ | 100 | 100 | | 102 | 6,8* |
| *Blue River 36R19 | | 94 | 95 | | | 22* | *Croplan Genetics 4801VT3 | CB,CR,RR-GQJ | 102 | 102 | | 102 | 6,9 |
| *Blue River 41R00 | | 98 | 97 | | | 22* | *Croplan Genetics 5338VT3 | CB,CR,RR-GQJ | 104 | 104 | 102 | 102 | 6*,9*,18* |
| *Blue River 44R57 | | 99 | 102 | | | 22* | *Croplan Genetics 5891VT3 | CB,CR,RR-GQJ | 105 | 104 | 107 | 102 | 6,9,16* |
| *Blue River 46M96 | | 100 | | 103 | 96 | 18* | *Croplan Genetics 591VT3 | CB,CR,RR-GQJ | 108 | 107 | 106 | 102 | 7,13,16,18* |
| *Blue River 48B30 | | 102 | 102 | 104 | 96 | 18*,21* | *Croplan Genetics 6125VT3 | CB,CR,RR-GQJ | 109 | 109 | | 102 | 7* |
| *Blue River 53R57 | | 105 | 102 | | | 21* | *Croplan Genetics 645TS | CB,CR,RR-GNJ | 111 | | 110 | 102 | 14*,16* |
| Blue River 59L36 | | 109 | | 106 | 96 | 18 | *Croplan Genetics 6818VT3 | CB,CR,RR-GQJ | 114 | | 112 | 102 | 14 |
| *Carharts Blue Top CR1000H | CB,CR,LL,RR-LPDJ | 100 | 99 | 102 | 57 | 8,11,15,18* | *Croplan Genetics 6831RHXT | CB,CR,LL,RR-LPDJ | 112 | | 112 | 102 | 14*,16* |
| Carharts Blue Top CR105RR | RR-J | 105 | 103 | | 35 | 6,9 | *Croplan Genetics 7505VT3 | CB,CR,RR-GQJ | 115 | 113 | | 102 | 14* |
| Carharts Blue Top CR1185V | CB,CR,RR-GQJ | 85 | 87 | | 57 | 10,12 | *Croplan Genetics S4900CR | CR,RR-MQJ | 102 | | 104 | 102 | 15,18* |
| *Carharts Blue Top CR1857V | CB,CR,RR-GQJ | 90 | 90 | | 33 | 8,10,12* | *Croplan Genetics S6100VT3 | CB,CR,RR-MGQJ | 110 | | 108 | 102 | 14*,16*,18* |
| *Carharts Blue Top CR5000V | CB,CR,RR-GQJ | 95 | 95 | | 57 | 8*,11 | Dahlco Seeds 7810VT3 | CB,CR,RR-GQJ | 81 | 85 | | 16 | 10 |
| *Carharts Blue Top CR6103R | RR-J | 103 | | 102 | 57 | 15* | Dahlco Seeds 7871GTCBLL | CB,LL,RR-KDS | 87 | 88 | | 16 | 10 |
| Carharts Blue Top CR9680V | CB,CR,RR-GQJ | 98 | 96 | | 33 | 8,11 | *Dahlco Seeds 7881VT3 | CB,CR,RR-GQJ | 88 | 88 | | 16 | 10* |
| Channel 190-20R Brand | RR-J | 90 | 87 | | 115 | 10 | Dahlco Seeds 7931GTCBLL | CB,LL,RR-KDS | 93 | 94 | | 16 | 11 |
| Channel 193-45R Brand | RR-J | 93 | 96 | | 115 | 11 | Dahlco Seeds 7951GTCBLL | CB,LL,RR-KDS | 95 | 95 | | 16 | 11 |
| Channel 194-70VT3 Brand | CB,CR,RR-GQJ | 94 | 94 | | 115 | 11 | Dahlco Seeds 8001GTCBLL | CB,LL,RR-KDS | 100 | 98 | | 16 | 8 |
| Channel 195-46VT3 Brand | CB,CR,RR-GQJ | 95 | 96 | | 115 | 11 | Dahlco Seeds 8021VT3 | CB,CR,RR-GQJ | 102 | 103 | | 16 | 9 |
| Channel 197-14VT3 Brand | CB,CR,RR-GQJ | 97 | 98 | | 115 | 8 | *Dahlco Seeds 8041GTCBLL | CB,LL,RR-KDS | 104 | 104 | 101 | 16 | 9,15* |
| *Channel 199-55VT3 Brand | CB,CR,RR-GQJ | 99 | 98 | | 37 | 8* | Dahlco Seeds 8111VT3 | CB,CR,RR-GQJ | 111 | 106 | | 16 | 9 |
| Channel 200-22VT3 Brand | CB,CR,RR-GQJ | 100 | 98 | | 115 | 8 | Dahlco Seeds 8841VT3 | CB,CR,RR-GQJ | 84 | 87 | | 16 | 10 |
| *Channel 205-94VT3 Brand | CB,CR,RR-GQJ | 105 | | 103 | 37 | 18* | Dahlco Seeds 8961VT3 | CB,CR,RR-GQJ | 96 | 96 | | 16 | 11 |
| *Channel 207-07VT3 Brand | CB,CR,RR-GQJ | 107 | 109 | | 115 | 7* | *Dahlman D4356 | CB,LL,RR-KDH | 85 | 86 | | 100 | 12* |
| Channel 210-57VT3 Brand | CB,CR,RR-GQJ | 110 | 110 | | 115 | 7 | *Dahlman D4356CBLL | CB,LL-KD | 85 | 86 | | 100 | 10,12* |
| *Channel 210-61VT3 Brand | CB,CR,RR-GQJ | 110 | 110 | 109 | 37 | 7,16* | *Dahlman R4342VT3 | CB,CR,RR-GQJ | 86 | 87 | | 100 | 10*,12* |
| *Channel 211-84T Brand | CB,CR,RR-GNJ | 111 | | 111 | 37 | 14*,16 | Dahlman R4807VT3 | CB,CR,RR-GQJ | 95 | 95 | | 100 | 11 |
| *Channel 212-34VT3 Brand | CB,CR,RR-GQJ | 112 | | 113 | 37 | 14* | Dahlman R4928VT3 | CB,CR,RR-GQJ | 98 | 96 | | 100 | 11 |
| | | | | | | | *Dairyland HiDF3000-9 | CB,CR,RR-GQJ | 100 | | 102 | 33 | 15,18* |

† Code = Trait(Gene); B = bmr(bm3); C = IMI(IT); D = LL(T25); F, G, K, L = Bt-CB: (Bt176, Mon810, Bt11, TC1507); H, J, S = RR(MonGA21, Nk603, SYTGA21); M = Leafy(lfy1); N, P, Q, R = Bt-CR (Mon863, DAS591227, Mon88017, MIR604); X = Unknown.

Table 2. Companies and hybrids included in the 2009 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

| Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables | Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables |
|----------------------------|------------------|------------------------|-----|-----|------|------------|----------------------------|------------------|------------------------|-----|-----|------|------------|
| | | Co. | GRM | SRM | Trt. | | | | Co. | GRM | SRM | Trt. | |
| Dairyland HiDF3007 | | 107 | | 110 | 33 | 16 | Dyna Gro 57V98 | CB,CR,RR-GQJ | 110 | 109 | | 49 | 7 |
| Dairyland HiDF3008-4 | CR,RR-NJ | 108 | | 106 | 33 | 13 | Dyna Gro V2483VT3 | CB,CR,RR-GQJ | 84 | 86 | | 49 | 10,12 |
| * Dairyland HiDF3104 | | 104 | | 105 | 33 | 13*,15* | Dyna Gro V3593VT3 | CB,CR,RR-GQJ | 95 | 95 | | 49 | 8,11 |
| * Dairyland HiDF3105Q | CB,CR,LL,RR-LPDJ | 105 | | 108 | 33 | 13*,16* | Dyna Gro V3883VT3 | CB,CR,RR-GQJ | 98 | 98 | | 49 | 8 |
| * Dairyland HiDF3110-6 | RR-J | 110 | | 112 | 33 | 14* | * Dyna Gro V4993V3 | CB,CR,RR-GQJ | 109 | 107 | | 49 | 7* |
| * Dairyland HiDF3187-7 | CB,RR-GJ | 87 | | 90 | 33 | 19* | Foundation Direct 8900CBR | CB,CR-KR | 87 | 88 | | 102 | 10 |
| * Dairyland HiDF3195Q | CB,CR,LL,RR-LPDJ | 95 | | 93 | 37 | 17* | * Foundation Direct OR8484 | | 108 | 108 | | | 21* |
| * Dairyland ST6992 | RR-J | 92 | 94 | | 33 | 11* | Foundation Direct OR8520 | | 104 | 104 | | | 22 |
| Dairyland ST7286 | CB,RR-GJ | 86 | 87 | | 33 | 12 | * Foundation Direct OR8585 | | 108 | 102 | | | 21* |
| * Dairyland ST7789 | CB,RR-GJ | 89 | | 90 | 33 | 19* | * Foundation Direct OR8643 | | 102 | 103 | | | 21* |
| Dairyland ST7790 | CB,RR-LJ | 90 | 90 | | 37 | 12 | * Foundation Direct OR8772 | | 100 | 100 | | | 22* |
| * Dairyland ST8208 | CB,CR,LL-LPD | 108 | 105 | | 33 | 7* | Foundation Direct OR8777 | | 95 | 98 | | | 22 |
| Dairyland ST8404 | CB,CR,LL-LPD | 101 | | 101 | 33 | 15 | * Foundation Direct OR8797 | | 95 | 96 | | | 22* |
| Dairyland ST9006 | CB,CR,RR-GQJ | 106 | 105 | | 33 | 9 | * Foundation Direct OR8800 | | 94 | 96 | | | 22* |
| * Dairyland ST9009 | CB,CR,RR-GQJ | 111 | 109 | | 33 | 7* | Foundation Direct OR8870 | | 92 | 94 | | | 23 |
| Dairyland ST9196 | CB,CR,RR-GQJ | 96 | 95 | | 33 | 11 | * Foundation Direct OR8900 | | 87 | 90 | | | 23* |
| * Dairyland ST9206Q | CB,CR,LL,RR-LPDJ | 106 | 106 | | 37 | 7*,9 | * Foundation Direct OR8934 | | 85 | 84 | | | 23* |
| * Dairyland ST9395 | CB,CR,RR-GQJ | 95 | 94 | | 33 | 11* | * Foundation Direct OR8947 | | 91 | 89 | | | 23* |
| Dairyland ST9500Q | CB,CR,LL,RR-LPDJ | 99 | 99 | | 37 | 8,11 | Foundation Direct OR8951 | | 86 | 93 | | | 23 |
| * Dairyland ST9597Q | CB,CR,LL,RR-LPDJ | 97 | 97 | | 37 | 8*,11 | Foundation Direct OR8999 | | 90 | 91 | | | 23 |
| Dairyland ST9703Q | CB,CR,LL,RR-LPDJ | 103 | 104 | | 37 | 6,9 | * G2 Genetics 1H-005AHXLL | CB,LL-LD | 105 | 104 | | 102 | 6* |
| * Dairyland ST9789 | CB,CR,RR-GQJ | 89 | 87 | | 33 | 10*,12* | * G2 Genetics 1H-005HXLL | CB,LL-LD | 105 | 104 | | 102 | 6* |
| Dairyland ST9799 | CB,CR,RR-GQJ | 99 | 98 | | 33 | 8,11 | * G2 Genetics 1X-911HXTLL | CB,CR,LL-LPD | 111 | 109 | 112 | 102 | 7,14* |
| * Dairyland ST9810 | CB,CR,RR-GQJ | 110 | 109 | | 33 | 7* | * G2 Genetics 5H-005RRHX | CB,LL,RR-LDJ | 105 | 106 | | 102 | 6*,9 |
| Dekalb DKC33-54(RR2) | RR-J | 83 | 83 | | 7 | 12 | * G2 Genetics 5H-007ARRHX | CB,LL,RR-LDJ | 107 | 108 | | 102 | 7* |
| * Dekalb DKC35-19(RR2YGCB) | CB,RR-GJ | 85 | 85 | 86 | 7 | 12,19* | * G2 Genetics 5H-007RRHX | CB,LL,RR-LDJ | 107 | 106 | | 102 | 7*,9 |
| Dekalb DKC36-34(VT3) | CB,CR,RR-GQJ | 86 | 86 | | 7 | 12 | G2 Genetics 5H-199RRHX | CB,LL,RR-LDJ | 99 | 101 | | 102 | 6 |
| * Dekalb DKC38-89(VT3) | CB,CR,RR-GQJ | 89 | 89 | 93 | 7 | 10*,12,17* | * G2 Genetics 5H-501RRHX | CB,LL,RR-LDJ | 100 | 99 | | 102 | 6,8*,11 |
| Dekalb DKC40-20(VT3) | CB,CR,RR-GQJ | 90 | 89 | | 7 | 10 | G2 Genetics 5H-506RRHX | CB,LL,RR-LDJ | 106 | 107 | | 102 | 7,9 |
| * Dekalb DKC41-60(VT3) | CB,CR,RR-GQJ | 91 | 95 | 93 | 7 | 11,17* | * G2 Genetics 5H-797RRHX | CB,LL,RR-LDJ | 97 | 96 | 94 | 102 | 8*,11*,17* |
| Dekalb DKC43-27(VT3) | CB,CR,RR-GQJ | 93 | 94 | | 7 | 11 | G2 Genetics 5H-884RRHX | CB,LL,RR-LDJ | 84 | 86 | | 102 | 10 |
| Dekalb DKC44-92(RR2) | RR-J | 94 | 95 | | 7 | 11 | * G2 Genetics 5H-905RRHX | CB,LL,RR-LDJ | 105 | 104 | | 102 | 6* |
| * Dekalb DKC45-79(VT3) | CB,CR,RR-GQJ | 95 | 96 | 96 | 7 | 8,11,17* | * G2 Genetics 5H-999RRHX | CB,LL,RR-LDJ | 99 | 99 | | 102 | 6,8*,11 |
| Dekalb DKC46-60(VT3) | CB,CR,RR-GQJ | 96 | 95 | | 7 | 8,11 | G2 Genetics 5X-000RRHXT | CB,CR,LL,RR-LPDJ | 100 | 101 | | 102 | 6*,11 |
| * Dekalb DKC48-37(VT3) | CB,CR,RR-GQJ | 98 | 95 | 93 | 7 | 8,11,17* | G2 Genetics 5X-594RRHXT | CB,CR,LL,RR-LPDJ | 94 | 97 | | 102 | 11 |
| Dekalb DKC50-44(VT3) | CB,CR,RR-GQJ | 100 | 100 | | 7 | 6,8 | * G2 Genetics 5X-614RRHXT | CB,CR,LL,RR-LPDJ | 114 | | 112 | 102 | 14* |
| * Dekalb DKC52-59(VT3) | CB,CR,RR-GQJ | 102 | 101 | 103 | 7 | 6*,9*,18* | G2 Genetics 5X-707RRHXT | CB,CR,LL,RR-LPDJ | 107 | 108 | | 102 | 7* |
| Dekalb DKC54-49(VT3) | CB,CR,RR-GQJ | 104 | | 102 | 7 | 15 | * G2 Genetics 5X-711RRHXT | CB,CR,LL,RR-LPDJ | 111 | 109 | 112 | 102 | 7*,14* |
| * Dekalb DKC55-07(VT3) | CB,CR,RR-GQJ | 105 | 103 | 106 | 7 | 6,13*,16* | G2 Genetics 5X-802RRHXT | CB,CR,LL,RR-LPDJ | 102 | 103 | | 102 | 6 |
| * Dekalb DKC55-24(VT3) | CB,CR,RR-GNJ | 105 | | 105 | 7 | 13* | * G2 Genetics 5X-905RRHXT | CB,CR,LL,RR-LPDJ | 105 | 105 | | 102 | 6*,9 |
| * Dekalb DKC57-50(VT3) | CB,CR,RR-GQJ | 107 | 109 | 106 | 7 | 7*,13,16* | * G2 Genetics 5X-909BRRHXT | CB,CR,LL,RR-LPDJ | 109 | 108 | | 102 | 7* |
| * Dekalb DKC57-66(VT3) | CB,CR,RR-GQJ | 107 | 108 | | 7 | 7* | * G2 Genetics 5X-909RRHXT | CB,CR,LL,RR-LPDJ | 109 | 106 | 107 | 102 | 7*,9,13 |
| * Dekalb DKC59-64(VT3) | CB,CR,RR-GQJ | 109 | 109 | 107 | 7 | 7*,13* | * G2 Genetics 5X-915RRHXT | CB,CR,LL,RR-LPDJ | 115 | | 112 | 102 | 14* |
| * Dekalb DKC61-69(VT3) | CB,CR,RR-GQJ | 109 | 110 | 107 | 7 | 7*,13* | Garst 83C55-3000GT | CB,CR,LL,RR-KRDS | 108 | | 109 | 102 | 13 |
| * Dyna Gro 52V01 | CB,CR,RR-GQJ | 87 | 86 | | 49 | 12* | * Garst 84A53GTCBLL | CB,LL,RR-KDS | 107 | | 107 | 102 | 13* |
| * Dyna Gro 54V78 | CB,CR,RR-GQJ | 96 | 95 | | 49 | 11* | Garst 84F34-3000GT | CB,CR,LL,RR-KRDS | 105 | | 108 | 102 | 13 |
| Dyna Gro 55R10 | CB,LL,RR-LDJ | 101 | | 104 | 49 | 18 | Garst 86G35-3000GT | CB,CR,LL,RR-KRDS | 100 | | 104 | 102 | 15 |
| * Dyna Gro 55V48 | CB,CR,RR-GQJ | 102 | | 102 | 49 | 18* | * Garst 86T82-3000GT | CB,CR,LL,RR-KRDS | 100 | | 102 | 102 | 15*,18 |
| Dyna Gro 56R29 | CB,CR,LL,RR-LPDJ | 106 | 105 | | 49 | 9 | Garst 87G94-3000GT | CB,CR,LL,RR-KRDS | 94 | | 97 | 102 | 17 |
| * Dyna Gro 57V40 | CB,CR,RR-GQJ | 111 | | 111 | 49 | 14*,16* | * Garst 88B38GTCBLL | CB,LL,RR-GDS | 90 | | 91 | 102 | 19* |

† Code = Trait(Gene); B = **bmr**(**bm3**); C = **IMI**(**IT**); D = **LL**(**T25**); F, G, K, L = **Bt-CB**: (Bt176, Mon810, Bt11, TC1507); H, J, S = **RR**(MonGA21, Nk603, SYTGA21); M = **Leafy**(**lfy1**); N, P, Q, R = **Bt-CR** (Mon863, DAS591227, Mon88017, MIR604); X = **Unknown**.

Table 2. Companies and hybrids included in the 2009 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

| Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables | Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables |
|-----------------------------|------------------|------------------------|-----|-----|------|----------|----------------------------|------------------|------------------------|-----|-----|------|-----------------|
| | | Co. | GRM | SRM | Trt. | | | | Co. | GRM | SRM | Trt. | |
| * Garst 89K64GTCBLL | CB,LL,RR-KDS | 83 | | 85 | 102 | 19* | Jung 7344VT3 | CB,CR,RR-GQJ | 91 | 96 | | 20 | 11 |
| Garst 89M60GTCBLL | CB,CR,LL,RR-KRDS | 82 | | 88 | 102 | 19 | Jung 7385VT3 | CB,CR,RR-GQJ | 94 | 96 | | 20 | 11 |
| * Gold Country 9608VT3 | CB,CR,RR-GQJ | 96 | 94 | | 7 | 11* | * Jung 7447VT3 | CB,CR,RR-GQJ | 95 | 96 | | 20 | 8*,11 |
| | | | | | | | Jung 7454VT3 | CB,CR,RR-GQJ | 98 | 98 | | 20 | 8 |
| * Golden Harvest H64673000G | CB,CR,RR-KRS | 87 | 85 | 87 | 102 | 12,19* | * Jung 7475VT3 | CB,CR,RR-GQJ | 100 | 98 | | 20 | 8* |
| Golden Harvest H82113000G | CB,CR,RR-KRS | 105 | 104 | | 102 | 9 | * Jung 7487VT3 | CB,CR,RR-GQJ | 99 | 97 | | 20 | 8* |
| Golden Harvest H88183000G | CB,CR,RR-KRS | 110 | 109 | | 102 | 7 | * Jung 7514VT3 | CB,CR,RR-GQJ | 102 | 103 | | 20 | 9* |
| * Great Lakes 3730G3VT3 | CB,CR,RR-GQJ | 87 | 87 | | 20 | 10* | Kaltenberg K2615LLGTBt11 | CB,LL,RR-KDS | 82 | 84 | | 49 | 12 |
| * Great Lakes 4041G3VT3 | CB,CR,RR-GQJ | 90 | 88 | | 20 | 10* | * Kaltenberg K2951LLBt11 | CB,LL-KD | 85 | 87 | | 49 | 12* |
| Great Lakes 4481G3VT3 | CB,CR,RR-GQJ | 94 | 95 | 96 | 20 | 11,17 | Kaltenberg K3023RRPlus | CB,CR,RR-GQJ | 86 | 86 | | 49 | 12 |
| Great Lakes 4689G3VT3 | CB,CR,RR-GQJ | 95 | 98 | | 20 | 8 | * Kaltenberg K3038LLBt11 | CB,LL-KD | 87 | 88 | | 49 | 12* |
| * Great Lakes 5306G3VT3 | CB,CR,RR-GQJ | 96 | 99 | | 20 | 8* | * Kaltenberg K3039LLGTBt11 | CB,LL,RR-KDS | 87 | 88 | | 49 | 12* |
| * Great Lakes 5335X2 | CB,LL-LD | 103 | | 101 | 20 | 15* | Kaltenberg K3625LLGTBt11 | CB,LL,RR-KDS | 93 | 90 | | 49 | 12 |
| Great Lakes 5450 | | 103 | 103 | | 20 | 9 | Kaltenberg K3843VT3 | CB,CR,RR-GQJ | 95 | 96 | 93 | 49 | 8,11,19 |
| Great Lakes 5783RR | RR-J | 107 | 109 | | 20 | 7 | Kaltenberg K4013VT3 | CB,CR,RR-GQJ | 95 | 95 | | 49 | 8,11 |
| * Great Lakes 5939G3VT3 | CB,CR,RR-GQJ | 109 | 110 | | 20 | 7* | Kaltenberg K4053VT3 | CB,CR,RR-GQJ | 97 | 97 | | 49 | 8 |
| | | | | | | | Kaltenberg K4149LLGT3 | CB,CR,LL,RR-KRDS | 98 | 100 | | 48 | 8 |
| Growmark FS34SV3 | CB,CR,RR-GQJ | 84 | 86 | | 75 | 10 | Kaltenberg K4521LLRRHXT | CB,CR,LL,RR-LPDJ | 100 | 100 | | 48 | 8 |
| * Growmark FS42SV3 | CB,CR,RR-GQJ | 92 | 94 | 94 | 75 | 11*,17* | * Kaltenberg K5332GT | RR-S | 104 | 104 | 104 | 48 | 9*,18* |
| * Growmark FS44SV3 | CB,CR,RR-GQJ | 94 | 95 | 95 | 75 | 11,17* | Kaltenberg K5355LLGTBt11 | CB,LL,RR-KDS | 104 | 106 | | 48 | 9 |
| Growmark FS49SV3 | CB,CR,RR-GQJ | 99 | 98 | | 75 | 11 | Kaltenberg K5588LLRRHXT | CB,CR,LL,RR-LPDJ | 105 | 105 | | 48 | 9 |
| * Growmark FS51SV3 | CB,CR,RR-GQJ | 101 | 102 | 103 | 75 | 9,15,18* | * Kaltenberg K6046LLHXT | CB,CR,LL,RR-LPDJ | 108 | 108 | | 48 | 7* |
| * Growmark FS53SV3 | CB,CR,RR-GQJ | 103 | 103 | 102 | 75 | 9,15* | * Kaltenberg K6645LLGT3 | CB,CR,LL,RR-KRDS | 110 | 109 | | 48 | 7* |
| * Growmark FS54SV3 | CB,CR,RR-GQJ | 104 | 103 | 100 | 75 | 9*,15* | * Kaltenberg K6663VT3 | CB,CR,RR-GQJ | 110 | 109 | | 48 | 7* |
| * Growmark FS55SV3 | CB,CR,RR-GQJ | 105 | 104 | 106 | 75 | 9,16* | * Kaltenberg K8097LF | M | 97 | | 94 | 48 | 19* |
| * Growmark FS56SV3 | CB,CR,RR-GQJ | 106 | 105 | 106 | 75 | 9,13,16* | * Kaltenberg K8610LF | M | 110 | | 105 | 48 | 18* |
| * Growmark FS58SV3 | CB,CR,RR-GQJ | 108 | 110 | | 75 | 7* | Kruger K1286RR | RR-J | 86 | 85 | | 37 | 10,12 |
| * Growmark FS60AV3 | CB,CR,RR-GQJ | 110 | 109 | | 75 | 7* | * Kruger K2090RRYGCB | CB,RR-GJ | 90 | 86 | 88 | 37 | 10,12*,19* |
| Growmark FS60K32NDS | CR,RR-NJ | 110 | | 111 | 75 | 14 | Kruger K2381RRYGCB | CB,RR-GJ | 81 | 84 | | 37 | 12 |
| Growmark FS60S44 | CB,CR,LL,RR-KRDS | 110 | 109 | 112 | 75 | 7,14 | * Kruger K6006VT3 | CB,CR,RR-GQJ | 106 | 106 | 106 | 37 | 7*,9*,13*,16* |
| * Growmark FS61BV3 | CB,CR,RR-GQJ | 111 | 110 | 113 | 75 | 7*,14* | Kruger K6007TS | CB,CR,RR-GNJ | 107 | 106 | | 37 | 7,9 |
| | | | | | | | * Kruger K6010VT3 | CB,CR,RR-GQJ | 110 | 110 | 112 | 37 | 7*,14* |
| * Hughes 3309 | CB,CR,LL-KRD | 103 | 103 | 104 | 37 | 6*,13 | * Kruger K6093VT3 | CB,CR,RR-GQJ | 93 | 97 | 91 | 37 | 8,11,19* |
| * Hughes 3639 | | 103 | 103 | 104 | 37 | 6,13* | * Kruger K6097VT3 | CB,CR,RR-GQJ | 97 | 97 | 96 | 37 | 8,11,17* |
| * Hughes 3928GT | RR-S | 104 | 103 | 106 | 37 | 6,13* | * Kruger K6102VT3 | CB,CR,RR-GQJ | 102 | 102 | 101 | 37 | 6,9*,15*,18* |
| Hughes 4592VT3 | CB,CR,RR-GQJ | 105 | 105 | 106 | 37 | 6,13 | * Kruger K6200VT3 | CB,CR,RR-GQJ | 100 | 98 | | 37 | 6,8*,11* |
| * Hughes 5594-3000GT | CB,CR,LL,RR-KRDS | 107 | 108 | 106 | 37 | 7*,13 | * Kruger K6208VT3 | CB,CR,RR-GQJ | 108 | 107 | 107 | 37 | 7*,13,16* |
| | | | | | | | * Kruger K6295VT3 | CB,CR,RR-GQJ | 95 | 96 | 92 | 37 | 8*,11,17*,19* |
| Johnson Seeds 4806VT3 | CB,CR,RR-GQJ | 80 | 84 | | 49 | 12 | * Kruger K6298VT3 | CB,CR,RR-GQJ | 98 | 98 | 97 | 37 | 6,8,11,15*,17 |
| Johnson Seeds 4846RR | RR-J | 84 | 86 | | 49 | 12 | Kruger K6378VT3 | CB,CR,RR-GQJ | 78 | 84 | | 37 | 12 |
| * Johnson Seeds 4852CBLL | CB,LL-KD | 85 | 86 | | 15 | 12* | * Kruger K6385VT3 | CB,CR,RR-GQJ | 85 | 84 | 87 | 37 | 10*,12,19 |
| * Johnson Seeds 4869VT3 | CB,CR,RR-GQJ | 86 | 87 | | 49 | 10*,12* | * Kruger K6388VT3 | CB,CR,RR-GQJ | 88 | 86 | 87 | 37 | 10,12,19* |
| Johnson Seeds 4874CBLLRW | CB,CR,LL-KRD | 87 | 88 | | 49 | 10,12 | Kruger K6401VT3 | CB,CR,RR-GQJ | 101 | 103 | 103 | 37 | 6,9,15,18 |
| * Johnson Seeds 4902CBLL | CB,LL-KD | 90 | 89 | | 15 | 10*,12* | Kruger K6408VT3 | CB,CR,RR-GQJ | 108 | 108 | | 37 | 7 |
| * Johnson Seeds 4925RRBT | CB,RR-GJ | 92 | 94 | | 15 | 11* | * Kruger K6410VT3 | CB,CR,RR-GQJ | 110 | 110 | | 37 | 7* |
| * Johnson Seeds 4966VT3 | CB,CR,RR-GQJ | 96 | 94 | | 49 | 11* | * Kruger K6411VT3 | CB,CR,RR-GQJ | 111 | 110 | 112 | 37 | 7*,14 |
| Johnson Seeds 4977-3000GT | CB,CR,LL,RR-KRDS | 97 | 98 | | 49 | 11 | * Kruger K6412VT3 | CB,CR,RR-GQJ | 112 | 110 | 112 | 37 | 7,14* |
| Johnson Seeds 4983GTCBLL | CB,LL,RR-KDS | 98 | 97 | | 15 | 11 | * Kruger K6490VT3 | CB,CR,RR-GQJ | 90 | 88 | 91 | 37 | 10,12,19* |
| | | | | | | | * Kruger K6499VT3 | CB,CR,RR-GQJ | 99 | 100 | 97 | 37 | 6,8*,11,15*,17* |
| Jung 7236VT3 | CB,CR,RR-GQJ | 84 | 84 | | 20 | 10 | * Kruger K6606VT3 | CB,CR,RR-GQJ | 106 | 106 | 107 | 37 | 7*,9,13,16* |
| * Jung 7288VT3 | CB,CR,RR-GQJ | 86 | 86 | | 20 | 10* | | | | | | | |

† Code = Trait(Gene); B = **bmr**(bmr3); C = **IMI**(IT); D = **LL**(T25); F, G, K, L = **Bt-CB**: (Bt176, Mon810, Bt11, TC1507); H, J, S = **RR**(MonGA21, Nk603, SYTGA21); M = **Leafy**(lfy1); N, P, Q, R = **Bt-CR** (Mon863, DAS591227, Mon88017, MIR604); X = **Unknown**.

Table 2. Companies and hybrids included in the 2009 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

| Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables | Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables |
|----------------------------|------------------|------------------------|-----|-----|------|---------------|----------------------------|-------------------|------------------------|-----|-----|------|-------------|
| | | Co. | GRM | SRM | Trt. | | | | Co. | GRM | SRM | Trt. | |
| *Kusmaul GL802 | RR-S | 100 | 99 | | 119 | 8* | Michael Fields HM6 | | 106 | 108 | | | 21 |
| Kusmaul GL805 | RR-S | 100 | 104 | | 119 | 8 | | | | | | | |
| Kusmaul GL890GT | RR-H | 90 | 89 | | 119 | 10 | *Midwest Genetics 77125T | CB,CR,RR-GNJ | 110 | | 109 | 37 | 14* |
| *Kusmaul K810 | | 110 | 107 | | 119 | 20* | Munson 14236LL | LL-D | 87 | 88 | | 94 | 10 |
| Kusmaul KX-8002HXT | CB,CR,LL,RR-LPDJ | 102 | 101 | | 119 | 6 | Munson 18656LLRR | LL,RR-DJ | 100 | 101 | | 94 | 11 |
| *Kusmaul KX-8008HXT | CB,CR,LL,RR-LPDJ | 108 | 108 | | 119 | 7* | *Munson 5033CBLL | CB,LL-KD | 90 | 87 | | 94 | 10* |
| Kusmaul SB788VT3 | CB,CR,RR-GQJ | 88 | 87 | | 119 | 10 | Munson 5841HXTRR | CB,CR,RR-LPJ | 97 | 99 | | 94 | 11 |
| LG Seeds LG2358VT3 | CB,CR,RR-GQJ | 85 | 85 | | 75 | 12 | Mycogen 2C302 | CB,LL,RR-LDJ | 90 | 91 | | 37 | 10,12 |
| *LG Seeds LG2362VT3 | CB,CR,RR-GQJ | 86 | 88 | | 75 | 10* | Mycogen 2D140 | CB,CR,RR-GQJ | 83 | 85 | | 37 | 10,12 |
| *LG Seeds LG2411VT3 | CB,CR,RR-GQJ | 92 | 94 | 94 | 75 | 11,17* | Mycogen 2D503 | | 100 | 98 | | 37 | 8 |
| LG Seeds LG2426VT3 | CB,CR,RR-GQJ | 94 | 95 | | 75 | 11 | *Mycogen 2E696 | CB,CR,RR-GQJ | 109 | 109 | | 37 | 7* |
| LG Seeds LG2466VT3 | CB,CR,RR-GQJ | 95 | 96 | | 75 | 11 | *Mycogen 2J337 | CB,CR,RR-GQJ | 92 | 91 | | 37 | 11*,12 |
| LG Seeds LG2478VT3 | CB,CR,RR-GQJ | 97 | 97 | | 75 | 11 | Mycogen 2J463 | CB,CR,RR-GQJ | 96 | 96 | | 37 | 8,11 |
| LG Seeds LG2496VT3 | CB,CR,RR-GQJ | 100 | 98 | | 75 | 8 | *Mycogen 2K662 | CB,CR,LL,RR-LPDJ | 108 | 108 | | 37 | 7* |
| *LG Seeds LG2498BT | CB-G | 98 | 96 | | 75 | 8* | Mycogen 2M495 | CB,CR,RR-GQJ | 99 | 99 | | 37 | 8,11 |
| LG Seeds LG2509CL | C | 103 | | 101 | 75 | 15 | *Mycogen 2P484 | CB,CR,RR-GQJ | 97 | 95 | | 37 | 8*,11 |
| *LG Seeds LG2510VT3 | CB,CR,RR-GQJ | 103 | 103 | | 75 | 9* | Mycogen 2T220 | CB,CR,RR-GQJ | 86 | 87 | | 37 | 12 |
| LG Seeds LG2514VT3 | CB,CR,RR-GQJ | 103 | 104 | | 75 | 9 | Mycogen 2W587 | CB,CR,LL,RR-LPDJ | 105 | 105 | | 37 | 9 |
| *LG Seeds LG2532VT3 | CB,CR,RR-GQJ | 106 | 107 | | 75 | 7* | *Mycogen 2Y547 | CB,CR,RR-GQJ | 103 | 102 | | 37 | 6*,9* |
| LG Seeds LG2548RR | RR-J | 110 | 109 | | 75 | 7 | Mycogen F28726 | CB,LL,RR-LDH | 109 | | 108 | 37 | 16,18 |
| *LG Seeds LG2552VT3 | CB,CR,RR-GQJ | 110 | | 113 | 75 | 14* | Mycogen F29309 | CB,LL,RR-LDH | 95 | | 96 | 37 | 19 |
| Legacy Seeds L2650RR | RR-J | 84 | 86 | | 59 | 12 | *Mycogen F2F489 | CB,LL,RR-BLDH | 99 | | 98 | 37 | 17* |
| Legacy Seeds L2820GTCBLL | CB,LL,RR-KDS | 88 | 88 | | 43 | 12 | *Mycogen F2F622 | CB,LL,RR-BLDH | 109 | 109 | 37 | | 16*,18 |
| *Legacy Seeds L2850VT3 | CB,CR,RR-GQJ | 88 | 86 | | 59 | 12* | Mycogen F2F700 | CB,CR,LL,RR-BLPDH | 113 | 112 | 37 | | 16 |
| *Legacy Seeds L3009RRBT | CB,RR-GJ | 92 | 94 | 90 | 59 | 11*,19* | Mycogen F2F797 | B | 115 | 113 | 37 | | 14 |
| *Legacy Seeds L3333VT3 | CB,CR,RR-GQJ | 95 | 94 | | 59 | 11* | Mycogen TMF2L418 | CB,CR,LL,RR-LPDH | 94 | | 96 | 37 | 17,19 |
| Legacy Seeds L3538VT3 | CB,CR,RR-GQJ | 95 | 97 | | 59 | 8,11 | *Mycogen TMF2Q296 | RR-H | 86 | | 91 | 37 | 19* |
| Legacy Seeds L4009HXTRR | CB,CR,LL,RR-LPDJ | 100 | 100 | 104 | 43 | 8,18 | *Mycogen TMF2Q716 | CB,CR,LL,RR-LPDH | 109 | 106 | 37 | | 16*,18* |
| *Legacy Seeds L4258VT3 | CB,CR,RR-GQJ | 102 | 104 | 102 | 59 | 9,18* | *Mycogen TMF2Q759 | CB,CR,LL,RR-LPDH | 113 | 112 | 37 | | 14* |
| Legacy Seeds L4938VT3 | CB,CR,RR-GQJ | 105 | 103 | | 59 | 6 | *Mycogen TMF2R521 | CB,CR,RR-GNH | 102 | 102 | 37 | | 15*,18* |
| Legacy Seeds L5309GT | RR-S | 104 | 104 | 104 | 43 | 6,15 | Mycogen TMF2W587 | CB,CR,LL,RR-LPDH | 104 | 103 | 37 | | 15 |
| *Legacy Seeds L5350GTCBLL | CB,LL,RR-KDH | 104 | 103 | 102 | 43 | 6,9*,15 | *Mycogen X29301 | CB,CR,LL,RR-LPDJ | 96 | 97 | | 37 | 8*,11 |
| *Legacy Seeds L6609HXTRR | CB,CR,LL,RR-LPDJ | 108 | | 109 | 43 | 16* | Mycogen X29421 | CB,CR,LL,RR-LPDJ | 99 | 99 | | 37 | 8,11 |
| *Legend Seeds LR9703VT3 | CB,CR,RR-GQJ | 103 | 102 | 101 | 52 | 9,11,15*,18* | Mycogen X29507 | CB,CR,LL,RR-LPDJ | 105 | 104 | | 37 | 6 |
| *Legend Seeds LR9707RRHXT | CB,CR,LL,RR-LPDJ | 107 | 104 | 105 | 52 | 9,11,16*,18* | *Mycogen X29583 | CB,CR,LL,RR-LPDJ | 103 | 104 | | 37 | 6,9* |
| *Legend Seeds LR9787GTCBLL | CB,LL,RR-KDS | 87 | 90 | 91 | 15 | 8,10,15*,17* | Next Generation Seed 3787 | CB,LL,RR-KDH | 87 | 88 | | 43 | 12 |
| *Legend Seeds LR9993VT3 | CB,CR,RR-GQJ | 93 | 95 | 94 | 52 | 8*,11,15*,17* | Next Generation Seed 4893 | CB,LL,RR-KDS | 93 | 94 | | 43 | 11 |
| *Legend Seeds LR9996VT3 | CB,CR,RR-GQJ | 96 | 96 | 98 | 52 | 8*,11,15,17 | Next Generation Seed 4998 | CB,LL,RR-KDS | 98 | 97 | | 43 | 11 |
| *Lemke 2020 | | 85 | 86 | | 56 | 12* | Next Generation Seed 4999 | CB,LL,RR-KDS | 99 | 98 | | 43 | 11 |
| *Lemke 3081Bt | CB-G | 90 | 89 | | 56 | 10* | *Next Generation Seed 5804 | CB,LL,RR-KDH | 104 | 103 | 103 | 43 | 9*,18* |
| Lemke 4097VT3 | CB,CR,RR-GQJ | 95 | 94 | | 56 | 11 | *Next Generation Seed 5808 | | 108 | 107 | 107 | 43 | 13*,16*,20* |
| Lemke 5077 | CB,CR,RR-GQJ | 100 | 99 | | 56 | 8 | Next Generation Seed 5902 | CB,LL,RR-KDS | 102 | 105 | | 43 | 9 |
| *Lemke 5141HX | CB,CR,LL,RR-LPDJ | 98 | 98 | | 56 | 8* | Next Generation Seed 5903 | CB,LL,RR-KDS | 103 | 103 | | 43 | 6 |
| *Lemke 6110 | | 105 | | 108 | 56 | 16* | NuTech 0C-603AYGCB | CB-G | 103 | 103 | | 6 | 6 |
| Lemke 6121 | LL-D | 105 | 105 | | 56 | 9 | NuTech 0C-603YGCB | CB-G | 103 | 102 | | 6 | 6 |
| *Masters Choice OG530 | | 105 | 108 | | | 21* | NuTech 1B-202CBLL | CB,LL-KD | 102 | 102 | | 102 | 6 |
| *Masters Choice SP590 | | 116 | | 112 | 117 | 14* | *NuTech 1B-290CBLL | CB,LL-KD | 90 | 87 | | 102 | 10* |
| | | | | | | | NuTech 1N-001CBLLRW | CB,CR,LL-KRD | 101 | 103 | | 6 | 6 |
| Michael Fields HM3 | | 105 | 104 | | | 21 | NuTech 1N-695CBLLRW | CB,CR,LL-KRD | 95 | 97 | | 102 | 11 |

† Code = Trait(Gene); B = **bmr**(bmr3); C = **IMI**(IT); D = **LL**(T25); F, G, K, L = **Bt-CB**: (Bt176, Mon810, Bt11, TC1507); H, J, S = **RR**(MonGA21, Nk603, SYTGA21); M = **Leafy**(lfy1); N, P, Q, R = **Bt-CR** (Mon863, DAS591227, Mon88017, MIR604); X = **Unknown**.

Table 2. Companies and hybrids included in the 2009 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

| Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables | Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables |
|-----------------------|------------------|------------------------|-----|-----|------|------------|----------------------------|------------------|------------------------|-----|-----|------|--------------|
| | | Co. | GRM | SRM | Trt. | | | | Co. | GRM | SRM | Trt. | |
| NuTech 1N-887CBLLRW | CB,CR,LL-KRD | 87 | 88 | | 6 | 10 | *OMG 4L15 | | 106 | 106 | | 2 | 20* |
| NuTech 3A-383RR | RR-J | 83 | 86 | | 6 | 10 | OMG 4L92 | | 108 | 108 | | 2 | 20 |
| *NuTech 3A-804GT | RR-H | 104 | 102 | 100 | 6 | 6,15* | *OMG 4M62 | | 104 | 106 | | 2 | 20* |
| NuTech 3C-115RRYGCB | CB,RR-GJ | 115 | | 113 | 102 | 14 | OMG 5L96 | | 109 | 109 | | 2 | 20 |
| NuTech 3C-389GTCB | CB,LL,RR-KDS | 88 | 87 | | 102 | 10 | *OMG 6E11 | | 109 | 107 | | 2 | 20* |
| *NuTech 3C-889RRYGCB | CB,RR-GJ | 89 | 88 | 91 | 102 | 10*,17*,19 | | | | | | | |
| NuTech 3P-191RRYGPL | CB,CR,RR-GQJ | 90 | 90 | | 6 | 10 | *Organic UTC | | 105 | 101 | | | 21*,22*,23* |
| NuTech 3P-494RRYGPL | CB,CR,RR-GQJ | 94 | 95 | | 6 | 11 | *Organic UTC-Hand Weed | | 105 | 102 | | | 21*,22*,23* |
| *NuTech 3T-013VT3 | CB,CR,RR-GQJ | 113 | | 112 | 102 | 14* | | | | | | | |
| NuTech 3T-098AVT3 | CB,CR,RR-GQJ | 98 | 98 | | 6 | 11 | *Pioneer 33F88 | CB,CR,LL,RR-LPDJ | 114 | | 111 | 98 | 14*,16 |
| *NuTech 3T-098VT3 | CB,CR,RR-GQJ | 98 | 98 | 96 | 6 | 11,17* | *Pioneer 33W84 | CB,CR,LL,RR-LPDJ | 111 | 108 | | 98 | 7* |
| *NuTech 3T-106VT3 | CB,CR,RR-GQJ | 106 | 108 | | 7 | 7* | *Pioneer 34A89 | CB,CR,LL,RR-LPDJ | 109 | | 106 | 98 | 13*,16* |
| NuTech 3T-110VT3 | CB,CR,RR-GQJ | 110 | 109 | | 6 | 7 | *Pioneer 35F40 | CB,LL,RR-LDJ | 105 | 102 | 102 | 98 | 6*,9,18 |
| NuTech 3T-295VT3 | CB,CR,RR-GQJ | 95 | 95 | | 7 | 11 | *Pioneer 35F44 | CB,CR,LL,RR-LPDJ | 105 | | 106 | 98 | 13*,16*,18 |
| NuTech 3T-300VT3 | CB,CR,RR-GQJ | 100 | 98 | | 6 | 11 | *Pioneer 35K04 | CB,CR,LL,RR-LPDJ | 106 | 106 | | 98 | 7*,9 |
| *NuTech 3T-390VT3 | CB,CR,RR-GQJ | 90 | 89 | | 102 | 10* | *Pioneer 36V53 | CB,LL,RR-LDJ | 102 | 103 | 102 | 98 | 6*,9*,15,18* |
| NuTech 3T-393VT3 | CB,CR,RR-GQJ | 93 | 94 | | 6 | 11 | Pioneer 36Y26 | CB,CR,LL,RR-LPDJ | 101 | 100 | 99 | 98 | 9,11,18,19 |
| NuTech 3T-408VT3 | CB,CR,RR-GQJ | 108 | 109 | | 6 | 7 | Pioneer 37K11 | CB,LL,RR-LDJ | 99 | 98 | | 98 | 8,11 |
| *NuTech 3T-413VT3 | CB,CR,RR-GQJ | 113 | | 112 | 102 | 14* | *Pioneer 37Y14 | CB,CR,LL,RR-LPDJ | 99 | 96 | | 98 | 8*,11* |
| *NuTech 3T-484VT3 | CB,CR,RR-GQJ | 84 | 88 | | 6 | 10* | *Pioneer 38M60 | CB,CR,LL,RR-LPDJ | 94 | 90 | 90 | 98 | 11,12,19* |
| *NuTech 3T-600VT3 | CB,CR,RR-GQJ | 100 | 99 | 102 | 7 | 6,11,15* | *Pioneer 38N88 | CB,LL,RR-LDJ | 92 | 90 | 89 | 98 | 11,12,19* |
| NuTech 3T-601VT3 | CB,CR,RR-GQJ | 101 | 102 | | 7 | 6 | *Pioneer 38P43 | CB,CR,LL,RR-LPDJ | 95 | 94 | 90 | 98 | 11,17*,19 |
| NuTech 3T-603VT3 | CB,CR,RR-GQJ | 103 | 104 | 103 | 6 | 6,9,15 | *Pioneer 39B23 | CB,LL,RR-LDJ | 89 | 86 | | 98 | 12* |
| *NuTech 3T-706VT3 | CB,CR,RR-GQJ | 106 | 107 | | 7 | 7* | Pioneer 39V07 | CB,LL,RR-LDJ | 80 | 84 | | 98 | 12 |
| NuTech 3T-801VT3 | CB,CR,RR-GQJ | 101 | 103 | | 7 | 6 | | | | | | | |
| *NuTech 3T-894VT3 | CB,CR,RR-GQJ | 94 | | 94 | 102 | 17* | *Prairie Hybrids 0371 | | 92 | 92 | | | 21*,22,23* |
| *NuTech 3T-904VT3 | CB,CR,RR-GQJ | 104 | 103 | | 7 | 6* | *Prairie Hybrids 3081 | | 104 | 105 | | | 21* |
| NuTech 3U-216VTRRLFY | CR,RR-MQJ | 116 | | 113 | 102 | 14 | *Producers Hybrids 6360 | C | | | 102 | 99 | 15* |
| *NuTech 5B-398GTCBLL | CB,LL,RR-KDS | 98 | 98 | 95 | 102 | 11,17* | *Producers Hybrids 7077VT3 | CB,CR,RR-GQJ | 110 | | 110 | 99 | 16* |
| NuTech 5B-593GTCBLL | CB,CR,LL,RR-KRDS | 93 | 94 | | 102 | 11 | *Producers Hybrids 7325VT3 | CB,CR,RR-GQJ | 113 | | 110 | 99 | 16* |
| *NuTech 5B-804GTCBLL | CB,LL,RR-KDS | 104 | 103 | | 6 | 6,9* | | | | | | | |
| *NuTech 5B-887GTCBLL | CB,LL,RR-KDS | 87 | | 90 | 102 | 19* | *Renk RK212CBLL | CB,LL-KD | 82 | 85 | 84 | 99 | 12*,19* |
| NuTech 5N-809GTCBLLRW | CB,CR,LL,RR-KRDS | 109 | 108 | | 7 | 7,9 | *Renk RK268VT3 | CB,CR,RR-GQJ | 85 | 87 | 89 | 99 | 12,19* |
| *NuTech 5X-100RRHXT | CB,CR,LL,RR-LPDJ | 100 | | 103 | 102 | 18* | *Renk RK292GTCBLL | CB,LL,RR-KDS | 85 | 88 | 88 | 99 | 10,12,19* |
| NuTech 5X-805RRHXT | CB,CR,LL,RR-LPDJ | 105 | 105 | | 102 | 6 | *Renk RK302CBLL | CB,LL-KD | 89 | 89 | 92 | 99 | 10*,12,19* |
| | | | | | | | *Renk RK434RRYGCB | CB,RR-GJ | 92 | 95 | 91 | 99 | 8*,11*,19* |
| OBrien OB1019A | | 104 | 106 | | 102 | 20 | Renk RK501VT3 | CB,CR,RR-GQJ | 95 | 96 | 92 | 99 | 11,19 |
| OBrien OB1019S | | 104 | 106 | | 102 | 20 | *Renk RK563CBLLRW | CB,CR,LL-KRD | 97 | 98 | 95 | 99 | 8,11,17* |
| OBrien OB1040LLHXT | CB,CR,LL-KLPD | 104 | 105 | | 102 | 6 | *Renk RK570VT3 | CB,CR,RR-GQJ | 95 | 95 | | 99 | 8*,11 |
| OBrien OB1042LLHXT | CB,CR,LL-KLPD | 105 | 105 | | 102 | 6 | *Renk RK594GTCBLLRW | CB,LL,RR-KDS | 100 | 99 | 103 | 99 | 8*,18* |
| *OBrien OB1058 | | 105 | 105 | 106 | 102 | 13*,16*,20 | *Renk RK616VT3 | CB,CR,RR-GQJ | 99 | 96 | 96 | 99 | 8,11,17* |
| *OBrien OB1058S | | 105 | 105 | 106 | 102 | 13*,16*,20 | Renk RK670VT3 | CB,CR,RR-GQJ | 103 | 103 | | 99 | 6,9 |
| OBrien OB1087 | | 108 | 108 | | 102 | 20 | *Renk RK686VT3 | CB,CR,RR-GQJ | 104 | 103 | 102 | 99 | 6,9,15*,18* |
| *OBrien OB1090LLCB | CB,LL-KD | 109 | 110 | | 102 | 7* | *Renk RK698VT3 | CB,CR,RR-GQJ | 103 | 103 | 101 | 99 | 6,9,15*,18* |
| OBrien OB3089GT | RR-S | 89 | 88 | | 102 | 10 | Renk RK711RRHXTRA | CB,CR,LL,RR-LPDJ | 107 | 106 | 108 | 99 | 7,9,16 |
| OBrien OB859 | | 85 | 85 | | 102 | 10 | *Renk RK744VT3 | CB,CR,RR-GQJ | 107 | 107 | 107 | 99 | 7*,16* |
| *OBrien OB879LL | LL-D | 88 | 94 | | 102 | 8* | *Renk RK760VT3 | CB,CR,RR-GQJ | 106 | 106 | | 99 | 7* |
| OBrien OBX105HXTRR2 | CB,CR,LL,RR-LPDJ | 105 | 104 | | 102 | 6 | *Renk RK829VT3 | CB,CR,RR-GQJ | 111 | 109 | 112 | 99 | 7*,14* |
| OBrien OBX950LLCB | CB,LL-KD | 95 | 98 | | 102 | 8 | *Renk RK844VT3 | CB,CR,RR-GQJ | 110 | 109 | 110 | 99 | 7,14*,16* |
| | | | | | | | | | | | | | |
| OMG 3L19 | | 100 | 104 | | 2 | 20 | *Spangler LFT49 | M | 102 | | 102 | 99 | 13*,15* |
| OMG 4E97 | | 102 | 102 | | 2 | 20 | *Spangler LFT59 | M | 109 | | 106 | 99 | 13* |

† Code = Trait(Gene); B = **bmr**(bm3); C = **IMI**(IT); D = **LL**(T25); F, G, K, L = **Bt-CB**: (Bt176, Mon810, Bt11, TC1507); H, J, S = **RR**(MonGA21, Nk603, SYTGA21); M = **Leafy**(lfy1); N, P, Q, R = **Bt-CR** (Mon863, DAS591227, Mon88017, MIR604); X = **Unknown**.

Table 2. Companies and hybrids included in the 2009 trials. A star (*) indicates that the hybrid performed statistically similar to the highest hybrid for yield or performance index (P.I. or MILK2006) in one or more zones.

| Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables | Brand Hybrid | Traits - Genes | Relative Maturity Seed | | | | Tables |
|----------------|----------------|------------------------|-----|-----|------|------------|--------------------|------------------|------------------------|-----|-----|------|--------|
| | | Co. | GRM | SRM | Trt. | | | | Co. | GRM | SRM | Trt. | |
| Trelay 2RR530 | RR-J | 86 | 86 | 91 | 99 | 10,17 | *UW EX42 | | 112 | | 112 | 1 | 14* |
| Trelay 3T544 | CB,CR,RR-GQJ | 93 | 97 | | 99 | 11 | Unity 4490VT3 | CB,CR,RR-GQJ | 90 | 87 | | 7 | 10 |
| Trelay 4RR455 | RR-J | 95 | 94 | | 99 | 11 | Unity 4496VT3 | CB,CR,RR-GQJ | 96 | 94 | | 7 | 11 |
| *Trelay 4VT741 | CB,CR,RR-GQJ | 96 | 96 | | 99 | 8*,11 | Unity 4502VT3 | CB,CR,RR-GQJ | 102 | 106 | | 7 | 9 |
| *Trelay 5T128 | CB,CR,RR-GQJ | 101 | 103 | 102 | 99 | 9*,15*,18* | *Unity 4504VT3 | CB,CR,RR-GQJ | 104 | 103 | | 7 | 9* |
| *Trelay 5VT323 | CB,CR,RR-GQJ | 99 | 98 | | 99 | 8* | Viking O.5740 | | 104 | 102 | | | 22 |
| *Trelay 6T226 | CB,CR,RR-GQJ | 106 | 104 | 102 | 99 | 9,18* | *Viking O.6710 | | 98 | 97 | | | 22* |
| *Trelay 6VT981 | CB,CR,RR-GQJ | 107 | | 105 | 99 | 16*,18* | *Viking O.6901N | | 101 | 102 | | | 22* |
| Trelay 7T231 | CB,CR,RR-GQJ | 111 | | 112 | 99 | 14 | *Viking O.6999N | | 99 | 102 | | | 22* |
| *Trelay 8RR712 | RR-J | 117 | | 112 | 99 | 14* | *Viking O.7292 | | 96 | 95 | | | 22* |
| *UW EX27 | | 112 | | 112 | 1 | 14* | *Viking O.99-90N | | 90 | 97 | | | 21* |
| *UW EX31 | | 110 | | 109 | 1 | 16* | Wyffels W2681 | CB,CR,RR-GQJ | 101 | 103 | | 7 | 6 |
| *UW EX34 | | 110 | | 110 | 1 | 16* | Wyffels W3629 | CB,CR,LL,RR-KRDS | 104 | 103 | | 6 | 6 |
| *UW EX35 | | 110 | | 112 | 1 | 14* | Wyffels W5159 | CB,CR,LL,RR-KRDS | 108 | 109 | | 6 | 7 |
| *UW EX36 | | 110 | | 108 | 1 | 16* | *Wyffels W5281 | CB,CR,RR-GQJ | 108 | 107 | | 7 | 7* |
| *UW EX37 | | 110 | | 108 | 1 | 16* | Wyffels W6871 | CB,CR,RR-GQJ | 110 | 110 | | 7 | 7 |
| *UW EX38 | | 112 | | 111 | 1 | 14* | YIELDirect 4M07VT3 | CB,CR,RR-GQJ | 103 | 102 | | 6 | 6 |
| *UW EX39 | | 112 | | 112 | 1 | 14* | YIELDirect 4X34VT3 | CB,CR,RR-GQJ | 103 | 104 | | 7 | 6 |
| UW EX40 | | 112 | | 112 | 1 | 14 | | | | | | | |
| *UW EX41 | | 112 | | 112 | 1 | 14* | | | | | | | |

Table 3. List of seed treatments used on corn hybrids entered in the 2009 UW corn trials.

| Seed Trt. | Common Name | Seed Trt. | Common Name |
|-----------|-------------------------------------|-----------|------------------------------------|
| 1 | Captan | 56 | Cruiser+Maxim |
| 2 | Maxim | 57 | Maxim+Poncho250 |
| 6 | Cruiser | 59 | Poncho250+Trilex |
| 7 | Poncho250 | 75 | ApronXL+MaximXL+Poncho250+Trilex |
| 15 | Actellic+MaximXL | 94 | Allegiance+Captan+Trilex+Poncho250 |
| 16 | ApronXL+Maxim+Poncho1250 | 96 | Natural II |
| 20 | ApronFL+Maxim+Poncho250+Trilex | 98 | Dynasty+MaximXL+Poncho 250 |
| 33 | Allegiance+Captan+Poncho250 | 99 | MaximXL+Poncho250+Trilex |
| 35 | Allegiance+Captan+Lorsban+Poncho250 | 100 | Maxim+Poncho250+Trilex |
| 37 | ApronXL+Cruiser+Dynasty+MaximXL | 102 | Cruiser250 |
| 43 | Cruiser+MaximXL | 114 | ApronXL+Dynasty+MaximXL+Poncho250 |
| 48 | Actellic+Cruiser+MaximXL | 115 | ApronXL+Dynasty+MaximXL+Poncho1250 |
| 49 | Actellic+MaximXL+Poncho250 | 117 | Actellic+MaximXL+Poncho250+Trilex |
| 52 | Actellic+ApronXL+Maxim+Poncho250 | 119 | Captan+Cruiser250 |
| 54 | ApronXL+Maxim | | |

† Code = Trait(Gene); B = **bmr**(bm3); C = **IMI**(IT); D = **LL**(T25); F, G, K, L = **Bt-CB**: (Bt176, Mon810, Bt11, TC1507); H, J, S = **RR**(MonGA21, Nk603, SYTGA21); M = **Leafy**(lfy1); N, P, Q, R = **Bt-CR** (Mon863, DAS591227, Mon88017, MIR604); X = **Unknown**.